

MSstats.daily, version 2.1.6
Protein significance analysis for mass spectrometry-based
proteomics

Meena Choi, choi67@purdue.edu, www.msstats.org

March 28, 2014

Contents

| | | |
|----------|---|-----------|
| 1 | Statistical relative protein quantification: SRM, DDA and DIA experiments | 3 |
| 2 | Allowable data formats | 7 |
| 2.1 | SRM with stable isotope labeled reference peptides | 7 |
| 2.1.1 | 10-column format | 7 |
| 2.1.2 | Assigning the values of <code>Condition</code> , <code>BioReplicate</code> and <code>Run</code> | 8 |
| 2.1.3 | <code>MSnSet</code> format | 10 |
| 2.2 | Label-free DDA | 11 |
| 2.2.1 | 10-column format | 11 |
| 2.2.2 | Assigning the values of <code>Condition</code> , <code>BioReplicate</code> and <code>Run</code> | 11 |
| 2.2.3 | <code>MSnSet</code> format | 12 |
| 2.3 | Label-free DIA | 12 |
| 2.3.1 | 10-column format | 12 |
| 2.3.2 | Assigning the values of <code>Condition</code> , <code>BioReplicate</code> and <code>Run</code> | 12 |
| 2.3.3 | <code>MSnSet</code> format | 12 |
| 3 | Example: SRM with stable isotope labeled reference peptides, a time-course investigation of <i>S. Cerevisiae</i> | 13 |
| 3.1 | Experimental design | 13 |
| 3.2 | Reading the data | 13 |
| 3.3 | Pre-processing data and quality control of MS runs | 13 |
| 3.4 | Model-based inference | 18 |
| 3.4.1 | Setting up the model and testing proteins for differential abundance | 18 |
| 3.4.2 | Verifying the assumption of the model | 22 |
| 3.4.3 | Visualizing the results of protein-level tests for differential abundance | 24 |
| 3.5 | Sample size calculation for a future experiment | 29 |
| 3.5.1 | Minimal number of biological replicates per condition | 29 |
| 3.5.2 | Statistical power of a future experiment | 30 |
| 3.5.3 | Visualization of sample size calculations | 30 |
| 3.6 | Quantification of protein abundance in individual samples or conditions | 31 |

| | | |
|----------|---|-----------|
| 4 | Example: label-free DDA, a controlled spike-in experiment | 34 |
| 4.1 | Experimental design | 34 |
| 4.2 | Reading the data | 34 |
| 4.3 | Pre-processing data and quality control of MS runs | 34 |
| 4.4 | Model-based inference | 35 |
| 4.4.1 | Setting up the model and testing proteins for differential abundance . . | 35 |
| 4.4.2 | Verifying the assumption of the model | 37 |
| 4.4.3 | Visualizing the results of protein-level tests for differential abundance . | 37 |
| 4.5 | Sample size calculation for a future experiment | 38 |
| 4.6 | Quantification of protein abundance in individual samples or conditions | 38 |
| 5 | Example: label-free DIA, a group comparison study of <i>S.Pyogenes</i> | 39 |
| 5.1 | Experimental design | 39 |
| 5.2 | Reading the data | 39 |
| 5.3 | Pre-processing data and quality control of MS runs | 39 |
| 5.4 | Model-based inference | 40 |
| 5.4.1 | Setting up the model and testing proteins for differential abundance . . | 40 |
| 5.4.2 | Verifying the assumption of the model | 41 |
| 5.4.3 | Visualizing the results of protein-level tests for differential abundance . | 41 |
| 5.5 | Sample size calculation for a future experiment | 42 |
| 5.6 | Quantification of protein abundance in individual samples or conditions | 42 |
| 6 | Formal description of the statistical models | 43 |
| 6.1 | SRM with stable isotope labeled reference peptides | 43 |
| 6.2 | Label-free DDA | 46 |
| 6.3 | Label-free DIA | 46 |

1 Statistical relative protein quantification: SRM, DDA and DIA experiments

MSstats is an open-source R-based package for statistical relative quantification of peptides and proteins in mass spectrometry-based proteomic experiments. This document describes `MSstats.daily`, the most recent (development) version of the package, and its use through the command line.

Applicability

MSstats version 2.0 and above is applicable to multiple types of sample preparation, including label-free workflows, workflows that use stable isotope labeled reference proteins and peptides, and workflows that use fractionation. It is applicable to targeted Selected Reaction Monitoring (SRM), Data-Dependent Acquisition (DDA or shotgun), and Data-Independent Acquisition (DIA or SWATH-MS). It is applicable to experiments that make arbitrary complex comparisons of experimental conditions or times.

MSstats is currently not applicable to experiments that compare multiple metabolically labeled endogenous samples within a same run. It is not applicable to experiments with iTRAQ labeling. These experiments will be supported in the future.

Statistical functionalities

MSstats version 2.0 and above performs three analysis steps. The first step, *data processing and visualization*, transforms and normalizes the intensities of the peaks, and generates workflow-specific and customizable numeric summaries for data visualization and quality control.

The second step, *statistical modeling and inference*, automatically detects the experimental design (e.g. group comparison, paired design or time course, presence of labeled reference peptides or proteins) from the data. It then reflects the experimental design, the type of spectral acquisition strategy, and the scope of conclusions (e.g. restricted to the subjects, or expanded to the underlying populations), and fits an appropriate linear mixed model by means of `lm` and `lmer` functionalities in R. The model is used to detect differentially abundant proteins or peptides, or to summarize the protein or peptide abundance in a single biological replicate or condition (that can be used, e.g. as input to clustering or classification).

The third step, *statistical experimental design*, views the dataset being analyzed as a pilot study of a future experiment, utilizes the variance components of the current datasets, and calculates the minimal number of replicates necessary in the future experiment to achieve a pre-specified statistical power.

Interoperability with existing computational tools

MSstats takes as input data in a tabular .csv format, which can be generated by any spectral processing tool such as SuperHirn (Mueller *et al.*, 2007), MaxQuant (Cox and Mann, 2008), Progenesis, MultiQuant, OpenMS (Sturm *et al.*, 2008) or OpenSWATH.

For statistics experts, MSstats 2.0 and above satisfies the interoperability requirements of Bioconductor, and takes as input data in the MSnSet format (Gatto and Lilley, 2012). The

command line-based workflow is partitioned into a series of independent steps, that facilitate the development and testing of alternative statistical approaches. It complies with the maintenance and documentation requirements of Bioconductor.

Finally, **MSstats** 2.0 and above is available as an external tool within Skyline (MacLean *et al.*, 2010). The external tool support within Skyline manages **MSstats** installation, point-and-click execution, parameter collection in Windows forms and output display. Skyline manages the annotations of the experimental design, and the processing of raw data. It outputs a custom report, that is fed as a single stream input into **MSstats**. This design buffers proteomics users from the details of the R implementation, while enabling rigorous statistical modeling.

Availability

MSstats is available under the Artistic-2.0 license at msstats.org. **MSstats** as an external tool is available at <http://proteome.gs.washington.edu/software/Skyline/tools.html>. **MSstats** is now also available in Bioconductor (<http://www.bioconductor.org>). The development version of the package, called **MSstats.daily**, is most recent and is available at msstats.org. We suggest to use that if possible. The versioning of the main package is updated several times a year, to synchronise with the Bioconductor release.

Changes from previous versions implemented in **MSstats** 1.0 and **SRMstats**

For special cases of some experimental workflows, the underlying statistical methodology was previously described, and implemented in R-based packages **MSstats** 1.0 and **SRMstats** (Chang *et al.*, 2012; Clough *et al.*, 2012; Surinova *et al.*, 2013).

MSstats 2.0 supersedes **MSstats** 1.0 and **SRMstats**, in that it implements all the analysis steps that are available in these packages. In addition, it extends the methodology and the implementation, as follows.

- Unlike **MSstats** 1.0 (limited to label-free DDA experiments) and **SRMstats** (limited to SRM experiments), **MSstats** 2.0 integrates the statistical analysis steps across two sample preparation workflows (label-free, and workflow using labeled reference proteins or peptides), and three spectral acquisition strategies (SRM, DDA and DIA). The integration enables a greater flexibility of statistical modeling for each experiment type.
- **MSstats** includes new statistical capabilities:
 - *Data processing*: quantification of the extent of between-run interferences.
 - *Data visualization*: more flexible plots of the protein profiles using **ggplot2** functionalities in R, in particular displaying pre-specified proteins, customizing axis range, label and angle.
 - *Fitting linear mixed effects models*: fit appropriate linear models in specialized circumstances (e.g. proteins with a single replicate in a condition, proteins with a single feature, proteins with various patterns of missing intensities in groups or runs), custom imputation of missing values with low-intensity signals, and custom removal of spectral features. Fit appropriate models for experiments with and without technical replicates. Model unequal variability of spectral features using iterative weight least squares.

- *Diagnostics of the quality of model fit*: residual plots and Normal quantile-quantile plots across features of a protein, and separately for each feature of a protein to detect deviations from the model assumptions such as unequal variance.
 - *Calculation of the sample size*: support of multiple modeling options in the analysis of the future experiment, such as expanded or restricted scope of biological replicates, and experiments with or without systematic interferences.
 - *Summarization of protein abundance in a subject or in a condition on a relative scale*: support of label-free experiments and experiments with labeled reference proteins or peptides. Support of multiple output formats (long format and data matrix).
- **MSstats 2.0** and above facilitates the interoperability with existing computational tools. In addition to taking as input a table in a CSV format, it can now be used as an external tool with Skyline by researchers who are unfamiliar with R. It also supports input in the **MSnSet** format, and partitions the analysis into a series of separate well-defined steps for interoperability with Bioconductor.

Overview of the functionalities

The functionalities of **MSstats** are overviewed in Figure 1. To get started with **MSstats**, first load the package and then visit the help section of **MSstats**-package using the following code.

```
> library(MSstats.daily)
> ?MSstats.daily
```

Troubleshooting

To help troubleshoot potential problems with installation or functionalities of **MSstats**, a progress report is generated in a separate log file **msstats.log**. The file includes information on the R session (R version, loaded software libraries), options selected by the user, checks of successful completion of intermediate analysis steps, and warning messages. If the analysis produces an error, the file contains suggestions for possible reasons for the errors. If a file with this name already exists in working directory, a suffix with a number will be appended to the file name. In this way a record of all the analyses is kept. Please see the file “KnownIssues-Skyline-MSstatsV2.1.6.pdf” on the “Installation” page of msstats.org for a list of known issues and possible solutions for installation problem of **MSstats** external tool in Skyline.




Figure 1: Overview of the functionalities and of the associated functions in **MSstats** . Colored boxes indicate the actual function names.

2 Allowable data formats

2.1 SRM with stable isotope labeled reference peptides

2.1.1 10-column format

MSstats performs statistical analysis steps, that follow peak identification and quantitation. Therefore, input to MSstats is the output of other software tools (such as Skyline or MultiQuant) that read raw spectral files and identify and quantify spectral peaks. The preferred structure of data for use in MSstats is a .csv file in a “long” format with 10 columns. The file should contain the following variables: **ProteinName**, **PeptideSequence**, **PrecursorCharge**, **FragmentIon**, **ProductCharge**, **IsotopeLabelType**, **Condition**, **BioReplicate**, **Run**, **Intensity**. The variable names are fixed, but are not case-sensitive. This required input data is generated automatically when using MSstats report format in Skyline.

- (a) **ProteinName**: This column needs information about Protein id. Statistical analysis will be done separately for each unique label in this column. For peptide-level modeling and analysis, use peptide id in this column.
- (b)-(e) **PeptideSequence**, **PrecursorCharge**, **FragmentIon**, **ProductCharge**: The combination of these 4 columns defines a *feature* of a protein (in SRM experiments, it is a transition that is identified and quantified across runs). If the information for one or several of these columns is not available, please do not discard these columns but use a single fixed value across the entire dataset. For example, if the original raw data does not contain the information of **ProductCharge**, assign the value 0 to the entries in the column **ProductCharge** for the entire dataset. If the peptide sequences should be distinguished based on post-translational modifications, this column can be renamed to **PeptideModifiedSequence**. For example, this allows us to use the **PeptideModifiedSequence** column from the Skyline report.
- (f) **IsotopeLabelType**: This column indicates whether this measurement is based on the endogenous peptides (use “L”) or labeled reference peptides (use “H”).
- (g) **Condition**: For group comparison experiments, this column indicates groups of interest (such as “Disease” or “Control”). For time-course experiments, this column indicates time points (such as “T1”, “T2”, etc). If the experimental design contains both distinct groups of subjects and multiple time points per subject, this column should indicate a combination of these values (such as “Disease_T1”, “Disease_T2”, “Control_T1”, “Control_T2”, etc).
- (h) **BioReplicate**: This column should contain a unique identifier for each biological replicate in the experiment. For example, in a clinical proteomic investigation this should be a unique patient id. Patients from distinct groups should have distinct ids.

MSstats does not require the presence of technical replicates in the experiment. If the technical replicates are present, all samples or runs from a same biological replicate should have a same id. MSstats automatically detects the presence of technical replicates and accounts for them in the model-based analysis.

- (i) **Run:** This column contains the identifier of a mass spectrometry run. Each mass spectrometry run should have a unique identifier, regardless of the origin of the biological sample. In SRM experiments, if all the transitions of a biological or a technical replicate are split into multiple ‘methods’ due to the technical limitations, each method should have a separate identifier. When processed by Skyline, distinct values of runs correspond to distinct input file names. It is possible to use the actual input file names as values in the column **Run**.
- (j) **Intensity:** This column should contain the quantified signal of a feature in a run without any transformation (in particular, no logarithm transform). The signals can be quantified as the peak height or the peak of area under curve. Any other quantitative representation of abundance can also be used.

An example of an input dataset is shown in Figure 2. More details on assigning the values of **Condition**, **BioReplicate** and **Run**, depending on the structure of the experimental design, are given below.

| | A | B | C | D | E | F | G | H | I | J |
|----|-------------|-----------------|-----------------|-------------|---------------|------------------|-----------|--------------|------------|-----------|
| 1 | ProteinName | PeptideSequence | PrecursorCharge | FragmentIon | ProductCharge | IsotopeLabelType | Condition | BioReplicate | Run | Intensity |
| 2 | ACEA | EILGHEIFFDWELP | 3 y3 | | 0 H | | 1 ReplA | 1 | 66472.3847 | |
| 3 | ACEA | EILGHEIFFDWELP | 3 y3 | | 0 L | | 1 ReplA | 1 | 5764.16228 | |
| 4 | ACEA | EILGHEIFFDWELP | 3 y4 | | 0 H | | 1 ReplA | 1 | 101005.166 | |
| 5 | ACEA | EILGHEIFFDWELP | 3 y4 | | 0 L | | 1 ReplA | 1 | 61.65238 | |
| 6 | ACEA | EILGHEIFFDWELP | 3 y5 | | 0 H | | 1 ReplA | 1 | 90055.4993 | |
| 7 | ACEA | EILGHEIFFDWELP | 3 y5 | | 0 L | | 1 ReplA | 1 | 472.691803 | |
| 8 | ACEA | TDSEAATLISSTID | 2 y10 | | 0 H | | 1 ReplA | 1 | 43506.5425 | |
| 9 | ACEA | TDSEAATLISSTID | 2 y10 | | 0 L | | 1 ReplA | 1 | 217.203553 | |
| 10 | ACEA | TDSEAATLISSTID | 2 y7 | | 0 H | | 1 ReplA | 1 | 68023.0377 | |
| 11 | ACEA | TDSEAATLISSTID | 2 y7 | | 0 L | | 1 ReplA | 1 | 725.284308 | |
| 12 | ACEA | TDSEAATLISSTID | 2 y8 | | 0 H | | 1 ReplA | 1 | 68276.0489 | |
| 13 | ACEA | TDSEAATLISSTID | 2 y8 | | 0 L | | 1 ReplA | 1 | 243.658527 | |

Figure 2: Example dataset from an SRM experiment with stable isotope labeled reference peptides. The dataset is stored in a .csv file in a “long” format. Each row corresponds to a single intensity.

2.1.2 Assigning the values of Condition, BioReplicate and Run

The values of **Condition**, **BioReplicate**, **Run** depend on the design of the specific experiment.

1. Group comparison

In a group comparison design, the conditions (e.g., disease states) are profiled across non-overlapping sets of biological replicates (i.e. subjects). In this example there are 2 conditions, Disease and Control (in general the number of conditions can vary). There are 3 subjects (i.e. biological replicates) per condition (in general an equal number of replicates per condition is not required). Each subject has 2 technical replicate runs (in general technical replicates are not required, and their number per sample may vary). Overall, in this example there are $2 \times 3 \times 2 = 12$ mass spectrometry runs.

Table 1 shows the values of the columns **Condition**, **BioReplicate** and **Run** for this situation. It is important to note two things. First, the order of subjects and conditions in the experiment should be randomized, and run id does not need to represent the order of spectral acquisition. Second, the values of the columns are repeated for every

quantified transition. For example, if in each run the experiment quantifies 50 endogenous transitions and 50 labeled reference counterparts, then the input file has $12 \times 50 \times 2 = 1200$ lines. When a feature intensity is missing in a run, the data structure should contain a separate row for each missing value. The rows should include all the information (from **ProteinName** to **Run**), and indicate missing intensities with ‘NA’.

| Condition | BioReplicate | Run |
|-----------|--------------|-----|
| Disease | Subject1 | 1 |
| Disease | Subject1 | 2 |
| Disease | Subject2 | 3 |
| Disease | Subject2 | 4 |
| Disease | Subject3 | 5 |
| Disease | Subject3 | 6 |
| Control | Subject4 | 7 |
| Control | Subject4 | 8 |
| Control | Subject5 | 9 |
| Control | Subject5 | 10 |
| Control | Subject6 | 11 |
| Control | Subject6 | 12 |

Table 1: Possible values of the columns **Condition**, **BioReplicate** and **Run** in an experiment with a group comparison design.

2. Time course

The important feature of a time course experimental design is that a same subject (i. e. biological replicate) is repetitively measured across multiple time points. In this example there are 2 time points, Time1 and Time2 (in general the number of times can vary). There are 4 subjects (i.e. biological replicates) measured across times (in general an equal number of times per replicate is not required). There are no technical replicates (in general the number of technical replicates per sample may vary). Overall, in this example there are $2 \times 4 \times 1 = 8$ mass spectrometry runs.

Table 2 shows the values of the columns **Condition**, **BioReplicate** and **Run** for this situation. Comments on the order of the runs, on the number of lines in the input data structure, and on the handling of missing peak intensities are as in the group comparison design.

| Condition | BioReplicate | Run |
|-----------|--------------|-----|
| Time1 | Subject1 | 1 |
| Time2 | Subject1 | 2 |
| Time1 | Subject2 | 3 |
| Time2 | Subject2 | 4 |
| Time1 | Subject3 | 5 |
| Time2 | Subject3 | 6 |
| Time1 | Subject4 | 7 |
| Time2 | Subject4 | 8 |

Table 2: Possible values of the columns **Condition**, **BioReplicate** and **Run** in a time course experiment.

3. Paired design

Another frequently used experimental design is a *paired design*, where measurements from multiple conditions (such as healthy biopsy and disease biopsy) are taken from a same subject. The statistical model for this experimental design is the same as in the time course experiment, however the values in the columns of the input data may have a different appearance. In this example there are 2 subjects, PatientA and PatientB (in general the number of patients can vary). There are two conditions per subject, BiopsyHealthy and BiopsyTumor (in general the number of conditions per subject can exceed two). In this example there are 3 technical replicates of each type (in this example, the technical replicates are biopsies; in general these can also be replicate sample preparations or replicate mass spectrometry runs). Overall, in this example there are $2 \times 2 \times 3 = 12$ mass spectrometry runs.

Table 3 shows the values of the columns **Condition**, **BioReplicate** and **Run** for this situation. Comments on the order of the runs, on the number of lines in the input data structure, and on the handling of missing peak intensities are as in the group comparison design.

2.1.3 MSnSet format

MSstats also allows data to be in the format of **MSnSet**, which is consistent with the requirements of Bioconductor. The **MSnSet** format has several components, of which the most commonly accessed are **assayData**, **phenoData**, and **featureData**. **assayData** is a matrix of intensities, where each row corresponds a transition, and the columns correspond to sample ids. **phenoData** contains columns that describe the biological samples, conditions in the experiment. **featureData** contains columns describing the peptide features, such as the name or id of the underlying protein and information of features.

If the data are stored in the format **expressionSet**, group labels information is required. If more than one variable is listed in the argument **group**, then a concatenated variable is created based on all of the specified **group** variables. The remaining information (peptide feature ids, biological replicate ids, and abundance) can be extracted from the rows and columns of **featureData** and **phenoData**, or the users can assign them based on their experimental design.

| Condition | BioReplicate | Run |
|---------------|--------------|-----|
| BiopsyHealthy | PatientA | 1 |
| BiopsyHealthy | PatientA | 2 |
| BiopsyHealthy | PatientA | 3 |
| BiopsyTumor | PatientA | 4 |
| BiopsyTumor | PatientA | 5 |
| BiopsyTumor | PatientA | 6 |
| BiopsyHealthy | PatientB | 7 |
| BiopsyHealthy | PatientB | 8 |
| BiopsyHealthy | PatientB | 9 |
| BiopsyTumor | PatientB | 10 |
| BiopsyTumor | PatientB | 11 |
| BiopsyTumor | PatientB | 12 |

Table 3: Possible values of the columns **Condition**, **BioReplicate** and **Run** in an experiment with paired design.

2.2 Label-free DDA

2.2.1 10-column format

For label-free DDA experiments the required input is the 10-column format, the same as described in Section 2.1 for SRM experiments. In DDA experiments spectral features are defined as peptide ions, which are identified and quantified across runs. Since for label-free DDA experiments some of the columns **PeptideSequence**, **PrecursorCharge**, **FragmentIon**, and **ProductCharge** are not relevant, these columns will have a constant fixed value (such as “NA”) across the entire dataset. Furthermore, the column **IsotopeLabelType** will be set to “L” for the entire dataset. An example dataset is shown in Figure 3.

| ProteinName | PeptideSequence | PrecursorCharge | FragmentIon | ProductCharge | IsotopeLabelType | Condition | BioReplicate | Run | Intensity |
|-------------|-----------------|-----------------|-------------|---------------|------------------|-----------|--------------|-----|------------|
| bovine | S.PVDIDTK | 5 | NA | NA | L | C1 | 1 | 1 | 2636791.5 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C1 | 1 | 2 | 1992418.5 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C1 | 1 | 3 | 1982146.38 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C2 | 1 | 4 | 5019594 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C2 | 1 | 5 | 4560467.5 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C2 | 1 | 6 | 3627848.75 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C5 | 1 | 13 | 145511.83 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C5 | 1 | 14 | 291829.69 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C6 | 1 | 16 | 786667.38 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C6 | 1 | 17 | 705295.31 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C6 | 1 | 18 | 453448.78 |
| bovine | S.PVDIDTK | 5 | NA | NA | L | C3 | 1 | 7 | NA |

Figure 3: Example dataset from a label-free DDA experiment. The dataset is stored in a .csv file in a “long” format. Each row corresponds to a single intensity.

2.2.2 Assigning the values of Condition, BioReplicate and Run

Same as in Section 2.1.2.

2.2.3 MSnSet format

Same as in Section 2.1.3

2.3 Label-free DIA

2.3.1 10-column format

For label-free DIA experiments, the required input is the 10-column format, the same as described in Section 2.1 for SRM experiments. The values of the required columns can be extracted from the output of signal processing software such as Skyline or OpenSWATH. By default, the combination of the values in the columns **PeptideSequence**, **PrecursorCharge**, **FragmentIon**, **ProductCharge** uniquely identifies each spectral feature (i.e. a fragment ion identified and quantified across multiple runs). If the signal processing software does not provide the information on some of these columns but provides a unique feature identifier, it is possible to use this unique identifier instead of one of these columns. Furthermore, the column **IsotopeLabelType** is set to “L” for the entire dataset.

An example dataset is shown in Figure 4. In this example, feature id generated by OpenSWATH is used instead of **ProductCharge** to uniquely characterize each feature.

| ProteinName | PeptideSequence | PrecursorCharge | FragmentIon | ProductCharge | IsotopeLabelType | Condition | BioReplicate | Run | Intensity |
|-------------|------------------------------|-----------------|-------------|---------------|------------------|-----------|--------------|-----|-----------|
| 350748 | TPPAAVLLK | 2 | y7 | 109401 | L | 2 | 1 | 3 | 257486 |
| 350748 | TPPAAVLLK | 2 | y7 | 109401 | L | 2 | 2 | 4 | 141159 |
| 350748 | TPPAAVLLK | 2 | y7 | 109401 | L | 1 | 1 | 1 | 452908 |
| 350748 | TPPAAVLLK | 2 | y7 | 109401 | L | 1 | 2 | 2 | 348222 |
| 515084 | NIC[160]VNAIAPGFIESDMTGVLPEK | 3 | y3 | 7717 | L | 2 | 1 | 3 | 12753 |
| 515084 | NIC[160]VNAIAPGFIESDMTGVLPEK | 3 | y3 | 7717 | L | 2 | 2 | 4 | 12857 |
| 515084 | NIC[160]VNAIAPGFIESDMTGVLPEK | 3 | y3 | 7717 | L | 1 | 1 | 1 | 89652 |
| 515084 | NIC[160]VNAIAPGFIESDMTGVLPEK | 3 | y3 | 7717 | L | 1 | 2 | 2 | 76724 |
| 515084 | MVNEAIESLGSIDVLVNNAGITNDK | 3 | y9 | 57971 | L | 2 | 1 | 3 | 2052 |
| 515084 | MVNEAIESLGSIDVLVNNAGITNDK | 3 | y9 | 57971 | L | 2 | 2 | 4 | 1050 |
| 515084 | MVNEAIESLGSIDVLVNNAGITNDK | 3 | y9 | 57971 | L | 1 | 1 | 1 | 10772 |
| 515084 | MVNEAIESLGSIDVLVNNAGITNDK | 3 | y9 | 57971 | L | 1 | 2 | 2 | 10516 |

Figure 4: Example dataset from a label-free DIA experiment. The dataset is stored in a .csv file in a “long” format. Each row corresponds to a single intensity.

2.3.2 Assigning the values of Condition, BioReplicate and Run

Same as in Section 2.1.2.

2.3.3 MSnSet format

Same as in Section 2.1.3

3 Example: SRM with stable isotope labeled reference peptides, a time-course investigation of *S. Cerevisiae*

3.1 Experimental design

This is a subset of the dataset published by [Picotti *et al.* \(2009\)](#). The study targeted 45 proteins in the glycolysis/gluconeogenesis/TCA cycle/glyoxylate cycle network of *S. Cerevisiae*. Three biological replicates were analyzed at ten time points (T1-T10), while the organism transited through exponential growth in a glucose-rich medium (T1-T4), diauxic shift (T5-T6), post-diauxic phase (T7-T9), and stationary phase (T10). Prior to trypsinization, the samples were mixed with an equal amount of proteins from the same N15-labeled yeast sample, which was used as a reference. The goal was to detect changes in protein abundance across time points. Transcriptional activity under the same experimental conditions has also been previously investigated by (DeRisi *et al.*, 1997).

For reasons of space the full collection of 45 proteins is not included with the package, however it is available at msstats.org in a tabular format. Measurements for two of the targeted proteins are included as part of the package: protein IDHC (gene name IDP2) is differentially abundant between time points T1 and T7. Protein PMG2 (gene name GPM2) is not differentially abundant between these time points. The protein names are based on the SwissProt nomenclature. In `MSstats` this dataset is stored in a data structure `SRMRawData`. For details visit the help file using the following code.

```
> ?SRMRawData
```

3.2 Reading the data

The dataset is stored in the “long” format, as a data.frame labeled `SRMRawData` that can be accessed once the package is installed and loaded in R.

```
> head(SRMRawData)
```

| | ProteinName | PeptideSequence | PrecursorCharge | FragmentIon | ProductCharge |
|-----|-------------|-----------------|-----------------|-------------|---------------|
| 243 | IDHC | ATDVIVPEEGELR | 2 | y7 | NA |
| 244 | IDHC | ATDVIVPEEGELR | 2 | y7 | NA |
| 245 | IDHC | ATDVIVPEEGELR | 2 | y8 | NA |
| 246 | IDHC | ATDVIVPEEGELR | 2 | y8 | NA |
| 247 | IDHC | ATDVIVPEEGELR | 2 | y9 | NA |
| 248 | IDHC | ATDVIVPEEGELR | 2 | y9 | NA |

| | IsotopeLabelType | Condition | BioReplicate | Run | Intensity |
|-----|------------------|-----------|--------------|-----|-------------|
| 243 | H | 1 | ReplA | 1 | 84361.08350 |
| 244 | L | 1 | ReplA | 1 | 215.13526 |
| 245 | H | 1 | ReplA | 1 | 29778.10188 |
| 246 | L | 1 | ReplA | 1 | 98.02134 |
| 247 | H | 1 | ReplA | 1 | 17921.29255 |
| 248 | L | 1 | ReplA | 1 | 60.47029 |

3.3 Pre-processing data and quality control of MS runs

(1) Data processing steps and options

Possible data processing steps include:

- Logarithm transform with base 2 (default) or 10 of the intensities
- Normalization to remove systematic bias between mass spectrometry runs. The normalization is applied after the logarithm transform. For SRM experiments with stable isotope labeled reference peptides, the normalization is typically based on labeled reference peptides of all the proteins. There are several options for normalization.
 - With the option `normalization=FALSE`, no normalization is performed.
 - With the option `normalization="equalizeMedians"` (default), constant normalization shifts all the intensities in a run by a constant, to equalize the *median of reference intensities* across runs.
 - With the option `normalization="quantile"`, quantile normalization ([Amaratunga and Cabrera, 2001](#)) applies a non-linear transformation to all the intensities in a run, to equalize the *distribution of reference intensities* across runs.
 - With the option `normalization="globalStandardProtein"`, the normalization is applied to endogenous intensities. First, the normalization equalizes endogenous intensities of global standard proteins across runs. Second, it applies the same between-run shifts to the remaining endogenous proteins in the experiment. For this normalization, global standard proteins or peptides should be assigned in `nameStandards` option, such as `nameStandards=c("ABC","DEF")`.
- Calculation of between-run interference score of a feature. The score is defined as the correlation of the feature intensity across runs, and the mean intensity of the corresponding peptide across runs.
- Addition of incomplete rows in the input. `MSstats` requires that the input contains a separate row for every feature in every run. If `MSstats` detects incomplete rows, it will output the list of problematic features. The incomplete rows can be filled according to two options.
 - With the option `fillIncompleteRows=FALSE` (default), the error message will be reported and the data processing will be stopped. This is done to allow the users to check the data.
 - With the option `fillIncompleteRows=TRUE`, the same warning message will be shown, but then the incomplete rows will be filled in based on the best possible guess, while adding `intensity=NA`.

Overall, this step produces an output summarizing the experimental design. The warning messages are sent both to the console and to the log file, notifying the user. A warning message includes the list of problematic features, subjects, conditions and their labels (reference or endogenous). Therefore users can verify that `MSstats` guessed the nature of the incomplete rows correctly. In addition, `MSstats` can distinguish duplicate rows, which is multiple rows for a same feature in a same run which are sometimes generated by signal processing tools. In this case, the user needs to decide which rows should be

used. Moreover, new columns are added to the dataset, for use in the downstream statistical modeling and model-based inference. For example, variable **abundance** represents the intensity of the feature used in the statistical modeling. Depending on the signal processing options, the intensity may be on the log or on the normalized log scale.

If all the transitions in a biological or technical replicate are split into multiple ‘methods’ (and are recorded in multiple files), this structure of the data is detected automatically by MSstats, by reading the values of the column **Run**. In this case the normalization is performed separately for each method.

To get started, visit the help file using the following code.

```
> ?dataProcess
```

In the example dataset, the processing step is as follows.

```
> QuantData<-dataProcess(SRMRawData)
```

```
Summary of Features :
               count
# of Protein      2
# of Peptides/Protein  2-2
# of Transitions/Peptide 3-3

Summary of Samples :
               1 2 3 4 5 6 7 8 9 10
# of MS runs      3 3 3 3 3 3 3 3 3 3
# of Biological Replicates 3 3 3 3 3 3 3 3 3 3
# of Technical Replicates 1 1 1 1 1 1 1 1 1 1
```

```
> head(QuantData)
```

| | PROTEIN | PEPTIDE | TRANSITION | | FEATURE LABEL | | |
|----|----------------|---------------------|-----------------|----------------------------|---------------|----------------|-----------|
| 1 | IDHC | ATDVIVPEEGELR_2 | y7_NA | ATDVIVPEEGELR_2_y7_NA | H | | |
| 3 | IDHC | ATDVIVPEEGELR_2 | y8_NA | ATDVIVPEEGELR_2_y8_NA | H | | |
| 5 | IDHC | ATDVIVPEEGELR_2 | y9_NA | ATDVIVPEEGELR_2_y9_NA | H | | |
| 7 | IDHC | DQTNDQVTVDSATATLK_2 | y10_NA | DQTNDQVTVDSATATLK_2_y10_NA | H | | |
| 9 | IDHC | DQTNDQVTVDSATATLK_2 | y11_NA | DQTNDQVTVDSATATLK_2_y11_NA | H | | |
| 11 | IDHC | DQTNDQVTVDSATATLK_2 | y8_NA | DQTNDQVTVDSATATLK_2_y8_NA | H | | |
| | GROUP_ORIGINAL | SUBJECT_ORIGINAL | RUN | GROUP | SUBJECT | SUBJECT_NESTED | INTENSITY |
| 1 | 1 | ReplA | 1 | 0 | 0 | 0.0 | 84361.083 |
| 3 | 1 | ReplA | 1 | 0 | 0 | 0.0 | 29778.102 |
| 5 | 1 | ReplA | 1 | 0 | 0 | 0.0 | 17921.293 |
| 7 | 1 | ReplA | 1 | 0 | 0 | 0.0 | 4481.229 |
| 9 | 1 | ReplA | 1 | 0 | 0 | 0.0 | 1871.042 |
| 11 | 1 | ReplA | 1 | 0 | 0 | 0.0 | 2640.060 |
| | ABUNDANCE | METHOD | SuggestToFilter | | | | |
| 1 | 15.84764 | 1 | 0 | | | | |
| 3 | 14.34531 | 1 | 0 | | | | |
| 5 | 13.61273 | 1 | 0 | | | | |
| 7 | 11.61303 | 1 | 0 | | | | |

| | | | |
|----|----------|---|---|
| 9 | 10.35297 | 1 | 0 |
| 11 | 10.84970 | 1 | 0 |

(2) Visualization for explanatory data analysis

`dataProcessPlots` takes as input the quantitative data from the function `dataProcess`, and generates three types of plots for data visualization and quality control.

- *QC plot* (Figure 5) visualizes potential systematic biases between mass spectrometry runs. After constant normalization, the median intensities of reference transitions across all proteins should be equal between runs (Figure 5(b)). After quantile normalization, the distribution of reference intensities across all proteins should be equal between runs (Figure 5(c)). This step generates two types of QC plots: one for all the proteins combined, and the other separately for each protein (produced in a separate pdf file).
- *Profile plot* (Figure 6) helps identify potential sources of variation (both variation of interest and nuisance variation) for each protein. Such plots should be done after the normalization.
- *Condition plot* (Figure 7) visualizes potential systematic differences in protein intensities between conditions. With the option `scale=TRUE`, the levels of conditions are scaled according to their labels. If `scale=FALSE` (default), the conditions on the x-axis are equally spaced. With the option `interval="CI"` (default), error bars indicate the confidence interval with 0.95 significant level for each condition. With the option `interval="SD"`, error bars indicate the standard deviation for each condition. The intervals are for descriptive purposes only, as more refined model-based estimation is obtained as discussed below.

The plots have a number of layout options, including size and description of axes labels, output file name etc. For example, the option `which.Protein` can be used to make plots for a restricted subset of proteins of interest. The option `address` specifies the name of the folder storing pdf files with the plots. With the option `address=FALSE`, plots will be shown in the graphical window, but not saved in a file. If a file with this name already exists in working directory, a suffix with a number will be appended to the file name. In this way a record of all the analyses is kept.

To get started, visit the help file using the following code.

```
> ?dataProcessPlots
```

(a) Quality control plot

```
> dataProcessPlots(data=QuantData, type="QCPlot")
```





Figure 5: Quality control (QC) plots for all the proteins combined. X-axis: run. Y-axis: log-intensities of transitions. Reference/endogenous signals are in the left/right panel. (a) Before normalization (`normalization=FALSE`). (b) After constant normalization (`normalization="equalizeMedians"`). (c) After quantile normalization (`normalization="quantile"`). The plots visualize potential artifacts in mass spectrometry runs.

(b) Profile plot

```
> dataProcessPlots(data=QuantData,type="ProfilePlot")
```




Figure 6: Profile plots for proteins IDHC and PMG2 after normalization. X-axis: run. Y-axis: log-intensities of transitions. Reference/endogenous signals are in the left/right panel. Line colors indicate peptides and line types indicate transitions. The plots help identify potential sources of interesting and nuisance variation for each protein.

(c) Condition plot

```
> dataProcessPlots(data=QuantData,type="ConditionPlot")
```




Figure 7: Condition plots for Protein IDHC and PMG2. X-axis: condition. Y-axis: log ratio of endogenous over reference intensities of each transition in a run. Dots indicate the mean of log ratio for each condition. Error bars are confidence intervals with 0.95 significant level for each condition. The plots visualize the differences between conditions, which are of the main biological interest.

3.4 Model-based inference

3.4.1 Setting up the model and testing proteins for differential abundance

A statistical model formally characterizes the sources of variation of all the measurements that pertain to a protein, and helps distinguish the systematic patterns of differential abundance

from noise. It describes the relationship between a *response* (i.e. intensities of the observed transitions) and a set of variables that have been observed with the response. In proteomic experiments the variables may include spectral features, conditions under which replicates are observed (e.g., disease state, stress applied to the organism, or time point), and biological replicates. **MSstats** is based on a family of linear mixed-effects models. The formal description of the models is given in Section 6.

MSstats implements this statistical modeling in the function **groupComparison**. It supports several experimental designs, including *group comparison* where different biological replicates are measured in each condition, *time course*, where same biological replicates are measured at different time points, and *paired design*, where measurements with multiple conditions are acquired from each biological replicate. The experimental design is recognized automatically based on the structure of the input data. Other required modeling assumptions that should be specified by the users include:

1. *Labeling technique*: **labeled=TRUE** (default) reflects the presence of labeled reference peptides or proteins (specify **FALSE** for label-free experiments).
2. *Scope of biological replication*: **scopeOfBioReplication="expanded"** states that the model-based conclusions should be extended beyond the subjects selected for the study, and are at the level of the underlying biological populations. **scopeOfBioReplication="restricted"** (default) restricts the model-based conclusions to the subjects selected to the study.
3. *Scope of technical replication*: **scopeOfTechReplication="expanded"** (default) states that the model-based conclusions should be extended beyond the mass spectrometry runs performed in the study, and are at the level of the underlying populations of possible replicates of mass spectrometry runs. **scopeOfTechReplication="restricted"** restricts the model-based conclusions to the performed mass spectrometry runs. This latter option leads to more conservative tests for differential protein abundance. In the special case of balanced experimental designs (where the same number of intensities are observed for each combination of feature, subject and condition), tests with this option are equivalent to tests based on the log-ratio of intensities of endogenous and reference transitions.
4. *The interpretation of interferences in feature intensities*: **interference=TRUE** (default) indicates that the feature intensities are assumed to be subject to systematic interferences (such as post-translational modifications) that are reproducible across multiple repetitions of the experiment. **interference=FALSE** indicates that the quantified interferences are non-reproducible random artifacts that should be considered as noise.
5. *Unequal variance of spectral features*: **equalFeatureVar=TRUE** (default) states the assumption that all the features have equal noise variation between mass spectrometry runs. **equalFeatureVar=FALSE** states the assumption that different feature intensities have different associated variability. **MSstats** reflects this assumption, and fits this more flexible model using iteratively re-weighted least squares (Kutner *et al.* (2005)). With this procedure we: (a) fit the constant variance model, (b) fit a smooth relationship between the fitted expected value and the residual variance, and (c) re-fit the linear model using this smooth relationship as the weight. Lower intensities have lower weights in the

model-based conclusions (and intensities of a same feature can have different weights in different conditions).

6. *Handling excessive missing intensities:* When a feature of a protein is missing completely in a condition or in a MS run, parameters of the original statistical model may become unestimable for that protein, and an adjustment to the model is needed. In this case a warning message is sent to the console and to the log file. The user is notified of three possible actions: (1) `missing.action="nointeraction"` (default) indicates that the quantified interferences are non-reproducible random artifacts that should be considered as noise and has the same effect as `interference=FALSE`; (2) `missing.action="impute"` imputes the missing values indicated by "NA" with the average of the smallest normalized intensities across run, or (3) `missing.action="remove"` removes that feature from the dataset. The model refinements are only applied to the proteins with missing values.

The choice of the scope of biological and technical replication should be made prior to the analysis, based of the experimental goals. This choice leads to different statistical models, and different estimation procedures (based on `lmer` and `lm` in R). It generally yields different conclusions. For the remaining specifications, we recommend starting with `interference=TRUE` as it maximizes the sensitivity of the model-based conclusions, and `equalFeatureVar=TRUE` as it minimizes overfitting. These specifications can be subsequently refined based on the model-based diagnostics plots in Section 3.4.2. All the modeling choices are applied simultaneously to all the proteins in the experiment.

In addition to the modeling assumptions, `groupComparison` requires us to state the conditions that we would like to compare. The statistical model will be used to evaluate each protein for evidence of differential abundance between these conditions, while taking into account the experimental design, and the available sources of variation. The comparisons are specified using the option `contrast.matrix`. R stores the levels of conditions in the alphanumeric order, and therefore the command `levels(QuantData$GROUP_ORIGINAL)` is a useful tool for checking the names and the order of the conditions in the memory.

In the example dataset, suppose that we would like to test proteins for differential abundance between times T1 and T7. We test the null hypothesis *no change in protein abundance* against the alternative hypothesis *change in protein abundance*. In statistical terminology, this can be written as $H_0 : L = \mu_{T7} - \mu_{T1} = 0$ against the alternative $H_a : L = \mu_{T7} - \mu_{T1} \neq 0$, where μ_{T7} and μ_{T1} are the mean population abundances of the protein at times T7 and T1. Since the dataset has 10 time points, and since the time points are listed in the alphanumeric order, this comparison corresponds to the following coefficients associated with each time point: $(-1, 0, 0, 0, 0, 0, 1, 0, 0, 0)$. Taking these coefficients as input, `groupComparison` estimates the log fold change and the standard error of the difference in abundance, and performs the hypothesis testing separately for each protein. It then reports p-values adjusted for multiple testing across the entire protein set. The full code for this comparison is as follows.

```
> levels(QuantData$GROUP_ORIGINAL)

[1] "1" "2" "3" "4" "5" "6" "7" "8" "9" "10"

> comparison<-matrix(c(-1,0,0,0,0,0,1,0,0,0),nrow=1)
> row.names(comparison)<-"T7-T1"
```

In the line above, `row.names(comparison)` labels the comparison with a description (an arbitrary character string chosen by the user) that helps the readability of the results.

```
> testResultOneComparison<-groupComparison(contrast.matrix=comparison, data=QuantData)
> testResultOneComparison$ComparisonResult
```

| | Protein | Label | log2FC | SE | Tvalue | DF | pvalue | adj.pvalue |
|---|---------|-------|-----------|-----------|-----------|-----|-----------|------------|
| 1 | IDHC | T7-T1 | 6.0227709 | 0.1438273 | 41.875016 | 100 | 0.0000000 | 0.0000000 |
| 2 | PMG2 | T7-T1 | 0.3044394 | 0.2242354 | 1.357677 | 100 | 0.1776219 | 0.1776219 |

The result of the test for differential abundance is a table with columns **Protein**, **Label** (of the comparison), **log2 fold change (log2FC)**, **standard error of the log2 fold change (SE)**, **test statistic of the Student test (Tvalue)**, **degree of freedom of the Student test (DF)**, **raw p-values (pvalue)**, **p-values adjusted for comparisons across multiple proteins using the approach by Benjamini and Hochberg (adj.pvalue)**. The cutoff of the adjusted p-value corresponds to the cutoff of the False Discovery Rate ([Benjamini and Hochberg, 1955](#)). The positive values of **log2FC** indicate evidence in favor of $\mu_{T7} > \mu_{T1}$ (i.e. proteins upregulated in T7), while the negative values indicate evidence in favor of $\mu_{T7} < \mu_{T1}$ (i.e. proteins downregulated in T7), as compared to T1.

The same model can be used to perform several comparisons of conditions simultaneously. In the example dataset, suppose that we would like to test T3-T1, T7-T1, and T9-T1. This can be done as follows. The output of these steps is stored in **MSstats** in a data structure **testResultMultiComparisons**.

```
> comparison1<-matrix(c(-1,0,1,0,0,0,0,0,0,0),nrow=1)
> comparison2<-matrix(c(-1,0,0,0,0,0,0,1,0,0,0),nrow=1)
> comparison3<-matrix(c(-1,0,0,0,0,0,0,0,0,1,0),nrow=1)
> comparison<-rbind(comparison1,comparison2, comparison3)
> row.names(comparison)<-c("T3-T1","T7-T1","T9-T1")
> testResultMultiComparisons<-groupComparison(contrast.matrix=comparison,data=QuantData)
> testResultMultiComparisons$ComparisonResult
```

| | Protein | Label | log2FC | SE | Tvalue | DF | pvalue | adj.pvalue |
|---|---------|-------|------------|-----------|------------|-----|-----------|------------|
| 1 | IDHC | T3-T1 | 0.1052223 | 0.1438273 | 0.7315877 | 100 | 0.4661312 | 0.4661312 |
| 4 | PMG2 | T3-T1 | -0.1830632 | 0.2242354 | -0.8163883 | 100 | 0.4162186 | 0.4661312 |
| 2 | IDHC | T7-T1 | 6.0227709 | 0.1438273 | 41.8750159 | 100 | 0.0000000 | 0.0000000 |
| 5 | PMG2 | T7-T1 | 0.3044394 | 0.2242354 | 1.3576775 | 100 | 0.1776219 | 0.1776219 |
| 3 | IDHC | T9-T1 | 6.1204163 | 0.1438273 | 42.5539234 | 100 | 0.0000000 | 0.0000000 |
| 6 | PMG2 | T9-T1 | 0.0718434 | 0.2242354 | 0.3203927 | 100 | 0.7493392 | 0.7493392 |

In this example there are three comparisons. Therefore, **MSstats** adjusts the p-values separately among all the proteins in the first comparison, then among all the proteins in the second comparison, and then among all the proteins in the third comparison. To simultaneously control the overall FDR at the level, say, 0.1, set the FDR in each individual comparison to $0.1/3$.

The coefficients of the comparisons between conditions do not need to be integers. For example, the comparison between T1 and the average of T2 and T3 ($H_0 : T1 - \frac{1}{2}(T2 + T3)$) is expressed with the coefficients `c(1, -0.5, -0.5, 0, 0, 0, 0, 0, 0, 0)`.

To get started, visit the help file using the following code.

```
> ?groupComparison
```

3.4.2 Verifying the assumption of the model

Results based on the statistical models are accurate as long as the assumptions of the models hold. Here we focus on the assumption of the Normal distribution of the measurement errors, and also on the assumption of constant variance of the measurement errors (if this option is specified in the model above). The assumptions can be checked by examining the residuals of the model fit (i.e., the deviations of the observed intensities of the transition from their model-based predictions).

MSstats generates residual plots and Normal quantile-quantile plots such as in Figure 8. Figure 8(a) is an example of protein where the variance of the residuals is associated with the mean feature intensity. Figure 8(b) illustrates that such deviations from constant variance can be mistaken for deviations from Normality. The feature-specific quantile-quantile plots confirm that in this case the assumption of constant variance should be relaxed, however deviations from Normality are not a major concern.

In **MSstats** residual plots such as in Figure 8(a) can be produced for each protein in the dataset using the function `modelBasedQCPlots`, taking as input the results of model fitting and testing in Section 3.4.1, using option `type="ResidualPlots"`. An additional option `which.Protein` can be used to produce the plots for a subset of the proteins, and the option `address` can be used to store the resulting pdf file in a particular location. With the option `address=FALSE` the plots are shown in the graphical window.

The quantile-quantile plots can be produced for each protein in the dataset using the function `modelBasedQCPlots` with the option `type="QQPlots"`. With the option `feature.QQPlots="all"`, the quantile-quantile plot is produced with all features of a protein, as in Figure 8(b). With the option `feature.QQplots="byFeature"`, the quantile-quantile plots are produced separately by each feature of a proteins, also separately by reference and endogenous intensities, as in Figure 8(c) and (d). Only large deviations of transition intensities from the straight line are problematic.

To get started, visit the help file using the following code.

```
> ?modelBasedQCPlots
```

The code for producing these plots with the example dataset is as follows.

```
> modelBasedQCPlots(data=testResultMultiComparisons$ModelQC, type="ResidualPlots",  
+                   which.Protein="PMG2", address=FALSE)  
  
> modelBasedQCPlots(data=testResultMultiComparisons$ModelQC, type="QQPlots",  
+                   which.Protein="PMG2", feature.QQPlot="all", address=FALSE)  
  
> modelBasedQCPlots(data=testResultMultiComparisons$ModelQC, type="QQPlots",  
+                   which.Protein="PMG2", feature.QQPlot="byFeature", address=FALSE)
```




Figure 8: Plots for checking the assumption of constant variance of the measurement errors for protein PMG2. (a) Residual plot. X-axis: predicted log-intensity of the feature, on average over all the runs. Y-axis: observed minus predicted log-intensity. The features have unequal variance. The variance can also be viewed as function of mean intensity of the feature. (b)-(d): Normal quantile-quantile plots. X-axis: theoretical quantiles of the Normal(0,1) distribution. Y-axis: quantiles of the observed minus predicted log-intensities. Deviations from the straight line indicate deviations from the model assumptions. (b) All the features combined. (c) Separately for each feature, reference intensities only. (d) Separately for each feature, endogenous intensities only. Panels (c) and (d) indicate that for this protein the pattern is due to deviations from the assumption of constant variance, and not necessarily from the assumption of Normality. The features with lower intensities have a larger variance, and are more likely to deviate from the Normality assumption.

3.4.3 Visualizing the results of protein-level tests for differential abundance

The function `groupComparisonPlots` takes as input the results of model fitting and testing in function `groupComparison` in Section 3.4.1, and visualizes them with volcano plots, heatmaps, and comparison plots.

Volcano plots visualize the outcome of one comparison between conditions for all the proteins, and combines the information on statistical and practical significance. The y-axis displays the FDR-adjusted p-values on the negative \log_2 scale, and represents statistical significance. The horizontal dashed line represents the FDR cutoff. The points above the FDR cutoff line are statistically significant differentially abundant proteins. These points are colored in red for upregulated proteins, and in blue for downregulated proteins. The x-axis is the model-based estimate of log-fold change (the base of logarithm transform is the same as specified in the `logTrans` option of the `dataProcess` step), and represents practical significance. It is possible to specify a practical significance cutoff based on the estimate of fold change in addition to the statistical significance cutoff. If the fold change cutoff is specified, the points above the horizontal cutoff line but within the vertical cutoff line will be judged as not differentially abundant (and will be colored in black). The practical significance cutoff can only be applied in addition to the statistical significance cutoff (i.e. the fold change alone does not present enough evidence for differential abundance).

Figure 9 continues the example of comparing multiple time points in Section 3.4.1, and shows several representative volcano plots for the comparison T7-T1. Figure 9(a)-(b) summarizes the comparison T7-T1 for the proteins IDHC and PMG2 in the example dataset. Figure 9(a) is the default volcano plot, obtained while adjusting the size of the label font with `x.axis.size=18` and `y.axis.size=18`. Figure 9(b) illustrates the effect of specifying the fold-change `FCcutoff=70` and removing protein names. (This example is for illustration only, since the fold change cutoff is unrealistically high. In this example protein IDHC lost its significance status after applying the fold change cutoff.) The code for producing Figure 9(b) is as follows.

```
> groupComparisonPlots(data=testResultMultiComparisons$ComparisonResult,
+                       type="VolcanoPlot",FCcutoff=70,
+                       x.axis.size=18, y.axis.size=18,
+                       ProteinName=FALSE, which.Comparison=c("T7-T1"), address=FALSE)
```

Since volcano plots are most effective when the number of proteins is relatively large, Figure 9(c)-(d) shows similar plots for the entire collection of the proteins in the *S. Cerevisiae* investigation.

The plots have a number of layout options, including upper or lower limits of the axes, and output file name. The option `address` specifies the name of the folder storing pdf files with the plots. With the option `address=FALSE`, plots will be shown in the graphical window, but not saved in a file. If a file with this name already exists in working directory, a suffix with a number will be appended to the file name.




Figure 9: Volcano plot of the comparison T7-T1. X-axis: practical significance, model-based estimate of log-fold change. Y-axis: statistical significance, FDR-adjusted p-values on the negative log₂ scale. The dashed line represents the FDR cutoff (default `sig=0.05`). (a)-(b) Proteins IDHC and PMG2. (a) Default volcano plot (the size of the label font adjusted with `x.axis.size=18` and `y.axis.size=18`). (b) The effect of specifying the fold-change `FCcutoff=70` and removing protein names. (c)-(d) as in (a)-(b), but for all the 45 proteins in the experiment. In (d), `FCcutoff=1.5`.

Heatmaps illustrate the patterns of up- and down-regulation of proteins in several comparisons. Examples of heatmaps are shown in Figure 10. Columns in the heatmaps are comparisons of conditions, and rows are proteins. The heatmaps display signed FDR-adjusted p-values of the tests, colored in red/blue for significantly up-/down-regulated proteins, while taking into account the specified FDR cutoff and the additional optional fold change cutoff. Brighter colors indicate stronger evidence in favor of differential abundance. Black color represents proteins are not significantly differentially abundant.

The rows and columns of the heatmaps can be ordered with the option `clustering`, which performs hierarchical clustering with the Ward method (minimum variance). The option `clustering="protein"` (default) clusters the rows (proteins) in the space of comparisons, based on the values of $(\text{sign of comparison}) \cdot (-\log_2(\text{adjusted p-values}))$. The option `clustering="comparison"` clusters the columns in the space of proteins, based on the values of $(\text{sign of comparison}) \cdot (-\log_2(\text{adjusted p-value}))$. The option `clustering="both"` reorders both columns and rows.

The code for producing Figure 10(a)-(b) is as follows.


```
> groupComparisonPlots(data=testResultMultiComparisons$ComparisonResult,
+                       type="Heatmap", address=FALSE)

> groupComparisonPlots(data=testResultMultiComparisons$ComparisonResult,
+                       type="Heatmap", FCcutoff=70, address=FALSE)
```


Since heatmaps are most effective when the number of proteins is relatively large, Figure 10(c)-(d) show similar plots for the entire collection of the proteins in the *S. Cerevisiae* investigation.

When the number of proteins is very large, interpretations of the heatmap can become difficult. In this case it is possible to split the heatmap in multiple sub-figures. The option `numProtein` can be used to indicate the number of proteins plotted in a sub-figure.


The plots have a number of layout options, including size of axes labels, and output file name. The option `address` specifies the name of the folder storing pdf files with the plots. With the option `address=FALSE`, plots will be shown in the graphical window, but not saved in a file. If a file with this name already exists in working directory, a suffix with a number will be appended to the file name.




(a) 2 proteins in RawData, default options



(b) 2 proteins in RawData, FCcutoff=70



(c) 45 proteins, default options



(d) 45 proteins, FCcutoff=1.5




Figure 10: Heatmap of results of testing proteins for differential abundance in three pairwise comparisons of conditions. Columns: comparisons; rows: proteins. Red: significant up-regulation; blue: significant down-regulation; black: no significant change in abundance. Brighter colors indicate stronger evidence in favor of differential abundance. (a)-(b) Two proteins in the entire dataset. (a) FDR cutoff of 0.05, and no fold change cutoff. (b) FDR cutoff of 0.05 and fold change cutoff 70. (c)-(d) As in (a)-(b), but for all the 45 proteins in the experiment. In (d), the fold change cutoff is **FCcutoff=1.5**.

Comparison plots illustrate model-based estimates of log-fold changes, and the associated uncertainty, in several comparisons of conditions for one protein. Two comparison plots for the two proteins in the example dataset are shown in Figure 11. X-axis is the comparison of interest. Y-axis is the log fold change. The dots are the model-based estimates of log-fold change, and the error bars are the model-based 95% confidence intervals (the option `sig` can be used to change the significance level of significance). For simplicity, the confidence intervals are adjusted for multiple comparisons within protein only, using the Bonferroni approach. For proteins with N comparisons, the individual confidence intervals are at the level of $1 - \text{sig}/N$.

The plots have a number of layout options, including size of axes labels, lower and upper limits of the axes, output file name etc. For example, the option `which.Protein` can be used to make plots for a restricted subset of proteins of interest. The option `address` specifies the name of the folder storing pdf files with the plots. With the option `address=FALSE`, plots will be shown in the graphical window, but not saved in a file. If a file with this name already exists in working directory, a suffix with a number will be appended to the file name. In this way a record of all the analyses is kept.

```
> groupComparisonPlots(data=testResultMultiComparisons$ComparisonResult,
+                       type="ComparisonPlot")
```

To get started, visit the help file using the following code.

```
> ?groupComparisonPlots
```




Figure 11: Comparison plots illustrate model-based estimates of log-fold changes, and the associated uncertainty, in comparisons of conditions for a single protein. X-axis: comparisons of interest. Y-axis: log-fold change. Dots: model-based estimates of log-fold change. Error bars: model-based 95% confidence intervals. Horizontal line: log fold change=0 (i.e. no significant change).

3.5 Sample size calculation for a future experiment

This last analysis step views the dataset as a pilot study of a future experiment, utilizes its variance components, and calculates the minimal number of replicates required in a future experiment to achieve the desired statistical power.

The calculation is performed by the function `designSampleSize`, which takes as input the same data structure `QuantData` as the statistical modeling step in Section 3.4.1. The function fits the same statistical model to the data as in Section 3.4.1 and has the same options. The only exception is that the option `scopeOfTechReplication` is fixed to `scopeOfTechReplication="restricted"` to provide conservative results.

Sample size calculation assumes same experimental design (i.e. group comparison, time course or paired design) as in the current dataset, and uses the model fit to estimate the median variance components across all the proteins. Finally, sample size calculation assumes that a large proportion of proteins (specifically, 99%) will not change in abundance in the future experiment. This assumption also provides conservative results.

Using the estimated variance components, the function relates the number of biological replicates per condition (`numSample`, rounded to 0 decimal), average number of peptides per protein (`numPep`), average number of transitions per peptide (`numTran`), average statistical power across all the proteins (`power`), minimal fold change that we would like to detect (can be specified as a range, e.g. `desiredFC=c(1.1, 2)`), and the False Discovery Rate (FDR). The user should specify all these quantities but one, and the function will solve for the remainder. The quantity to solve for should be set to `= TRUE`.

To get started, visit the help section of `designSampleSize` using the following code.

```
> ?designSampleSize
```

3.5.1 Minimal number of biological replicates per condition

The code and the output below is an example for calculating the minimal number of biological replicates for a future experiment, while using as a pilot study the two proteins in the time-course investigation of *S. Cerevisiae*. In this example we solve for the number of biological replicates per condition (`numSample=TRUE`), while considering a range of the desired fold changes between 1.25 and 1.75 (`desiredFC=c(1.25,1.75)`).

```
> designSampleSize(data=QuantData,numSample=TRUE,numPep=3,numTran=4,power=0.8,
+                   desiredFC=c(1.25,1.75),FDR=0.05)
```

| | desiredFC | numSample | numPep | numTran | FDR | power | CV |
|----|-----------|-----------|--------|---------|------|-------|-------|
| 1 | 1.250 | 15 | 3 | 4 | 0.05 | 0.8 | 0.004 |
| 2 | 1.275 | 12 | 3 | 4 | 0.05 | 0.8 | 0.005 |
| 3 | 1.300 | 11 | 3 | 4 | 0.05 | 0.8 | 0.006 |
| 4 | 1.325 | 9 | 3 | 4 | 0.05 | 0.8 | 0.007 |
| 5 | 1.350 | 8 | 3 | 4 | 0.05 | 0.8 | 0.007 |
| 6 | 1.375 | 7 | 3 | 4 | 0.05 | 0.8 | 0.008 |
| 7 | 1.400 | 6 | 3 | 4 | 0.05 | 0.8 | 0.009 |
| 8 | 1.425 | 6 | 3 | 4 | 0.05 | 0.8 | 0.009 |
| 9 | 1.450 | 5 | 3 | 4 | 0.05 | 0.8 | 0.011 |
| 10 | 1.475 | 5 | 3 | 4 | 0.05 | 0.8 | 0.011 |
| 11 | 1.500 | 4 | 3 | 4 | 0.05 | 0.8 | 0.013 |

| | | | | | | | |
|----|-------|---|---|---|------|-----|-------|
| 12 | 1.525 | 4 | 3 | 4 | 0.05 | 0.8 | 0.013 |
| 13 | 1.550 | 4 | 3 | 4 | 0.05 | 0.8 | 0.013 |
| 14 | 1.575 | 4 | 3 | 4 | 0.05 | 0.8 | 0.013 |
| 15 | 1.600 | 3 | 3 | 4 | 0.05 | 0.8 | 0.017 |
| 16 | 1.625 | 3 | 3 | 4 | 0.05 | 0.8 | 0.016 |
| 17 | 1.650 | 3 | 3 | 4 | 0.05 | 0.8 | 0.016 |
| 18 | 1.675 | 3 | 3 | 4 | 0.05 | 0.8 | 0.016 |
| 19 | 1.700 | 3 | 3 | 4 | 0.05 | 0.8 | 0.016 |
| 20 | 1.725 | 2 | 3 | 4 | 0.05 | 0.8 | 0.023 |
| 21 | 1.750 | 2 | 3 | 4 | 0.05 | 0.8 | 0.023 |

3.5.2 Statistical power of a future experiment

The code and the output below is an example for calculating the statistical power of a future experiment, while using as a pilot study the two proteins in the time-course investigation of *S. Cerevisiae*. In this example we solve for the statistical power (`power=TRUE`), while considering a range of the desired fold changes between 1.25 and 1.75 (`desiredFC=c(1.25,1.75)`).

```
> designSampleSize(data=QuantData,numSample=2,numPep=3,numTran=4,power=TRUE,
+                   desiredFC=c(1.25,1.75),FDR=0.05)
```

| | desiredFC | numSample | numPep | numTran | FDR | power | CV |
|----|-----------|-----------|--------|---------|------|-------|-------|
| 1 | 1.250 | 2 | 3 | 4 | 0.05 | 0.01 | 0.016 |
| 2 | 1.275 | 2 | 3 | 4 | 0.05 | 0.01 | 0.016 |
| 3 | 1.300 | 2 | 3 | 4 | 0.05 | 0.01 | 0.015 |
| 4 | 1.325 | 2 | 3 | 4 | 0.05 | 0.01 | 0.015 |
| 5 | 1.350 | 2 | 3 | 4 | 0.05 | 0.01 | 0.015 |
| 6 | 1.375 | 2 | 3 | 4 | 0.05 | 0.02 | 0.014 |
| 7 | 1.400 | 2 | 3 | 4 | 0.05 | 0.03 | 0.014 |
| 8 | 1.425 | 2 | 3 | 4 | 0.05 | 0.05 | 0.014 |
| 9 | 1.450 | 2 | 3 | 4 | 0.05 | 0.08 | 0.014 |
| 10 | 1.475 | 2 | 3 | 4 | 0.05 | 0.12 | 0.013 |
| 11 | 1.500 | 2 | 3 | 4 | 0.05 | 0.16 | 0.013 |
| 12 | 1.525 | 2 | 3 | 4 | 0.05 | 0.21 | 0.013 |
| 13 | 1.550 | 2 | 3 | 4 | 0.05 | 0.26 | 0.013 |
| 14 | 1.575 | 2 | 3 | 4 | 0.05 | 0.32 | 0.013 |
| 15 | 1.600 | 2 | 3 | 4 | 0.05 | 0.38 | 0.012 |
| 16 | 1.625 | 2 | 3 | 4 | 0.05 | 0.43 | 0.012 |
| 17 | 1.650 | 2 | 3 | 4 | 0.05 | 0.49 | 0.012 |
| 18 | 1.675 | 2 | 3 | 4 | 0.05 | 0.54 | 0.012 |
| 19 | 1.700 | 2 | 3 | 4 | 0.05 | 0.59 | 0.012 |
| 20 | 1.725 | 2 | 3 | 4 | 0.05 | 0.64 | 0.012 |
| 21 | 1.750 | 2 | 3 | 4 | 0.05 | 0.69 | 0.011 |


3.5.3 Visualization of sample size calculations

The calculated relationship between the number of biological replicates per condition (`numSample`), average number of peptides per protein (`numPep`), average number of transitions per peptide (`numTran`), average statistical power across all the proteins (`power`), minimal fold change that we would like to detect (`desiredFC`), and the False Discovery Rate (`FDR`) can be


visualized using the function `designSampleSizePlots`. The function takes as input the output of `designSampleSize`.

Two examples visualizing the minimal number of biological replicates per condition and the statistical power are given below.

```
> result.sample<-designSampleSize(data=QuantData, numSample=TRUE, numPep=3, numTran=4,
+                               power=0.8, desiredFC=c(1.25,1.75), FDR=0.05)
> designSampleSizePlots(data=result.sample)
```



```
> result.power<-designSampleSize(data=QuantData, numSample=2, numPep=3, numTran=4,
+                               power=TRUE, desiredFC=c(1.25,1.75), FDR=0.05)
> designSampleSizePlots(data=result.power)
```



To get started, visit the help section using the following code.

```
> ?designSampleSizePlots
```

3.6 Quantification of protein abundance in individual samples or conditions

Many downstream analysis steps (such as clustering or classification of individual samples in the space of their protein profiles) require summary values of protein abundance in each biological replicate or in each condition, on a relative scale that is comparable between runs.

MSstats enables such model-based summarization with the function **quantification**. The function takes as input the same data structure **QuantData** as the statistical modeling step in Section 3.4.1. This step fits the same linear mixed model as in **groupComparison**, with the difference that it specifies the reduced scope of biological replication for *sample quantification*.

The option **type="sample"** (default) performs *sample quantification*, i.e. it outputs the estimates of relative protein abundance in each biological replicate. In the special case of proteins quantified with a single transition **MSstats** gives a message that this estimate is not reliable. In presence of completely missing values in biological replicate, the estimates will be zero. The option **type="group"** performs *group quantification*, i.e. it outputs the estimates of relative protein abundance in each condition, summarized over the biological replicates. The summarization is performed separately for endogenous and labeled reference transitions. In presence of completely missing values in a condition, the estimates will be zero.

MSstats supports two output formats. The option **format="matrix"** (default) outputs an array where rows are proteins, and columns are conditions (for group quantification), or combinations of biological replicate and condition ids (for sample quantification). The option **format="long"** produces an array where each row corresponding to relative protein abundances, and columns are **Protein**, **Condition**, **LogIntensities** (and **BioReplicate** in the case of sample quantification).

The code and the output below is an example of quantification for the two proteins in the time-course investigation of *S. Cerevisiae*. Sample quantification outputs model-based estimates of protein abundance in each biological replicate within each of the ten time points. Group quantification outputs model-based estimates of protein abundance in each of the ten time points, summarized over all the biological replicates.

```
> subQuant<-quantification(QuantData)
> head(subQuant)
```

| | ReplA_1 | ReplA_2 | ReplA_3 | ReplA_4 | ReplA_5 | ReplA_6 | ReplA_7 | ReplA_8 |
|------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|----------|
| IDHC | 6.482258 | 7.144148 | 6.753797 | 6.891645 | 7.51246 | 9.955367 | 12.76964 | 12.94651 |
| PMG2 | 10.392707 | 9.959031 | 10.078554 | 10.520689 | 10.27148 | 9.909567 | 10.39021 | 10.14813 |
| | ReplA_9 | ReplA_10 | ReplB_1 | ReplB_2 | ReplB_3 | ReplB_4 | ReplB_5 | ReplB_6 |
| IDHC | 12.81644 | 12.94724 | 6.959746 | 6.84807 | 6.586199 | 6.86521 | 7.096384 | 9.86911 |
| PMG2 | 10.49330 | 10.26594 | 10.400018 | 10.72289 | 10.514320 | 10.23663 | 10.388234 | 10.38962 |
| | ReplB_7 | ReplB_8 | ReplB_9 | ReplB_10 | ReplC_1 | ReplC_2 | ReplC_3 | |
| IDHC | 12.77337 | 12.95800 | 12.88595 | 12.95984 | 6.916612 | 7.065317 | 7.154622 | |
| PMG2 | 10.55574 | 10.03953 | 10.28494 | 10.13845 | 10.433715 | 10.598824 | 10.034071 | |
| | ReplC_4 | ReplC_5 | ReplC_6 | ReplC_7 | ReplC_8 | ReplC_9 | ReplC_10 | |
| IDHC | 7.411565 | 6.554255 | 9.884292 | 12.74810 | 12.96304 | 12.82765 | 12.90613 | |
| PMG2 | 10.779834 | 10.287874 | 10.334494 | 11.15592 | 10.11160 | 10.21011 | 10.17977 | |
| | Ref | | | | | | | |
| IDHC | 12.884638 | | | | | | | |
| PMG2 | 9.107043 | | | | | | | |

```
> groupQuant<-quantification(QuantData, type="Group")
> head(groupQuant)
```

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
|------|----------|----------|---------|---------|----------|----------|----------|
| IDHC | 6.786205 | 7.019178 | 6.83154 | 7.05614 | 7.054366 | 9.902923 | 12.76371 |


```
PMG2 10.408813 10.426913 10.20898 10.51238 10.315861 10.211228 10.70063
      8          9          10
IDHC 12.95585 12.84335 12.93774
PMG2 10.09976 10.32945 10.19472
```

To get started, visit the help section using the following code.

```
> ?quantification
```

4 Example: label-free DDA, a controlled spike-in experiment

4.1 Experimental design

Mueller *et al.* (2007) described a controlled spike-in experiment, where 6 proteins, (horse myoglobin, bovine carbonic anhydrase, horse Cytochrome C, chicken lysozyme, yeast alcohol dehydrogenase, rabbit aldolase A) were spiked into a complex background in known concentrations in a latin square design. The experiment contained 6 mixtures, and each mixture was analyzed in label-free LC-MS mode with 3 technical replicates (resulting in the total of 18 runs). Each protein was represented by 7-21 peptides, and each peptide was represented by 1-5 transition. The raw data can be accessed from (http://prottools.ethz.ch/muellelu/web/Latin_Square_Data.php).

4.2 Reading the data

The dataset is stored as part of the package in the “long” format, as a data.frame labeled `DDARawData`. It can be accessed once the package is installed and loaded in R.

```
> head(DDARawData)
```

| | ProteinName | PeptideSequence | PrecursorCharge | FragmentIon | ProductCharge |
|---|------------------|-----------------|-----------------|-------------|---------------|
| 1 | bovine | S.PVDIDTK_5 | NA | 5 | NA |
| 2 | bovine | S.PVDIDTK_5 | NA | 5 | NA |
| 3 | bovine | S.PVDIDTK_5 | NA | 5 | NA |
| 4 | bovine | S.PVDIDTK_5 | NA | 5 | NA |
| 5 | bovine | S.PVDIDTK_5 | NA | 5 | NA |
| 6 | bovine | S.PVDIDTK_5 | NA | 5 | NA |
| | IsotopeLabelType | Condition | BioReplicate | Run | Intensity |
| 1 | L | C1 | 1 | 1 | 2636792 |
| 2 | L | C1 | 1 | 2 | 1992418 |
| 3 | L | C1 | 1 | 3 | 1982146 |
| 4 | L | C2 | 1 | 4 | 5019594 |
| 5 | L | C2 | 1 | 5 | 4560468 |
| 6 | L | C2 | 1 | 6 | 3627849 |

4.3 Pre-processing data and quality control of MS runs

(1) Data processing steps and options

The data processing steps are as in Section 3.3. However, label-free DDA experiments require adjustments in normalization steps, as compared to SRM experiments with isotope labeled reference peptides. The normalization options for label-free DDA experiments are as follows:

- The normalization is applied after the logarithm transform.
 - With the option `normalization=FALSE`, no normalization is performed.
 - With the option `normalization="equalizeMedians"` (default), constant normalization shifts all the endogenous intensities in a run by a constant, to equalize *the median endogenous intensities* across runs.

- With the option `normalization="quantile"`, quantile normalization ([Amaratunga and Cabrera, 2001](#)) applies a non-linear transformation to all the endogenous intensities in a run, to equalize the *distribution of the endogenous intensities* across runs.
- With the option `normalization="globalStandardProtein"`, the normalization equalizes the endogenous intensities of global standard proteins across runs, and then applies the same between-run shifts to the remaining endogenous features in the experiment. For this normalization, the ids of global standard proteins or peptides should be indicated in option `nameStandards`, e.g. `nameStandards=c("Protein1","Protein2")`.

The code and the output for the controlled spike-in experiment are given below.

```
> DDAquant<-dataProcess(DDARawData)

Summary of Features :
               count
# of Protein           6
# of Peptides/Protein 11-32
# of Transitions/Peptide 1-1

Summary of Samples :
               C1 C2 C3 C4 C5 C6
# of MS runs           3  3  3  3  3  3
# of Biological Replicates 1  1  1  1  1  1
# of Technical Replicates  3  3  3  3  3  3
```

(2) Visualization for explanatory data analysis

As in Section 3.3. The example code generating profile plots, quality control plots and condition plots are given below.

```
> dataProcessPlots(data=DDAquant,type="ProfilePlot", ylimDown=10, ylimUp=32,
+                 featureName="NA",x.axis.size=18,text.size=7,width=7,height=7,address=FALSE)
> dataProcessPlots(data=DDAquant,type="QCPlot",ylimDown=10, ylimUp=32,
+                 width=7, height=7,address=FALSE)
> dataProcessPlots(data=DDAquant,type="ConditionPlot",width=7,height=7,address=FALSE)
```

Figure 12 shows the example profile plot for one of the spiked proteins.

4.4 Model-based inference

4.4.1 Setting up the model and testing proteins for differential abundance

The steps of setting up the statistical model for label-free DDA experiments are as in Section 3.4.1. However, since these experiments are characterized by non-systematic interferences, unequal variances of the measurements, and a large number missing intensities (especially among low-intensity features), we recommend starting with different options. To express the fact that the interferences are non-systematic, use `interference=FALSE`. To express the fact that the variances are unequal between the features, use `equalFeatureVar=FALSE`.




Figure 12: Profile plots for a spiked protein (Lysozyme) in the controlled spike-in DDA experiment, after normalization. X-axis: run. Y-axis: log- intensities of endogenous features. Line colors indicate peptide ions. The hallmarks of label-free DDA experiments are non-systematic interferences, unequal variances of the measurements, and a large number missing intensities (especially among low-intensity features). The plot helps identify potential sources of interesting and nuisance variation for each protein.

Since this is a label-free experiment, use `labeled=FALSE`. Currently the option `equalFeatureVar=FALSE` is only available with the reduced scope of biological and technical replication.

In addition to specifying the probability model, the user should specify the comparisons of interest. This step is the same as in Section 3.4.1. Example of the code and of the output for the controlled spike-in experiment are given below. In this example we are interested in 6 pairwise comparisons of neighboring conditions (condition 2 vs condition 1, condition 3 vs condition 2, etc).

```
> comparison1<-matrix(c(-1,1,0,0,0,0),nrow=1)
> comparison2<-matrix(c(0,-1,1,0,0,0),nrow=1)
> comparison3<-matrix(c(0,0,-1,1,0,0),nrow=1)
> comparison4<-matrix(c(0,0,0,-1,1,0),nrow=1)
> comparison5<-matrix(c(0,0,0,0,-1,1),nrow=1)
> comparison6<-matrix(c(1,0,0,0,0,-1),nrow=1)
> comparison<-rbind(comparison1,comparison2,comparison3,comparison4,comparison5,comparison6)
> row.names(comparison)<-c("C2-C1","C3-C2","C4-C3","C5-C4","C6-C5","C1-C6")

> testResultDDA<-groupComparison(contrast.matrix=comparison,data=DDAquant,labeled=FALSE,
+                               interference=FALSE,equalFeatureVar=FALSE)
> head(testResultDDA$ComparisonResult)
```

| | Protein Label | log2FC | SE | Tvalue | DF | pvalue |
|----|-----------------|------------|-----------|------------|-----|--------------|
| 1 | bovine C2-C1 | 0.4630392 | 0.1218444 | 3.800251 | 102 | 2.461202e-04 |
| 7 | chicken C2-C1 | 0.8424595 | 0.2562116 | 3.288140 | 90 | 1.439771e-03 |
| 13 | cyc_horse C2-C1 | 1.1860546 | 0.1428324 | 8.303822 | 251 | 6.217249e-15 |
| 19 | myg_horse C2-C1 | -6.2774924 | 0.3167507 | -19.818403 | 87 | 0.000000e+00 |
| 25 | rabbit C2-C1 | 1.2803467 | 0.4238242 | 3.020938 | 218 | 2.820945e-03 |
| 31 | yeast C2-C1 | 0.6173339 | 0.2754126 | 2.241488 | 109 | 2.702025e-02 |

```

      adj.pvalue
1  4.922405e-04
7  2.159656e-03
13 1.865175e-14
19 0.000000e+00
25 3.385134e-03
31 2.702025e-02

```

4.4.2 Verifying the assumption of the model

As in Section 3.4.2. Example code for the controlled spike-in experiment is given below. An example of the resulting residual plot is in Figure 13.

```

> modelBasedQCPlots(data=testResultDDA$ModelQC,type="ResidualPlots",
+                   featureName=FALSE, address=FALSE)
> modelBasedQCPlots(data=testResultDDA$ModelQC,type="QQPlots",
+                   feature.QQPlot="byFeature",address=FALSE)

```




Figure 13: Residual plot for one spiked protein, Alcohol dehydrogenase. X-axis: predicted log-intensity of the feature, on average over all the runs. Y-axis: observed minus predicted log-intensity. This particular protein shows only a minor deviation from the constant variance assumption.

4.4.3 Visualizing the results of protein-level tests for differential abundance

As in Section 3.4.3. Example code for the controlled spike-in experiment is given below. An example of the resulting heatmap is in Figure 14.

```

> groupComparisonPlots(data=testResultDDA$ComparisonResult,
+                      type="VolcanoPlot",address=FALSE)
> groupComparisonPlots(data=testResultDDA$ComparisonResult,
+                      type="Heatmap",address=FALSE)
> groupComparisonPlots(data=testResultDDA$ComparisonResult,
+                      type="ComparisonPlot",address=FALSE)

```




Figure 14: Heatmap of results of 6 pairwise comparisons of conditions in the controlled spike-in experiment. Columns: comparisons; rows: proteins. Red: significant up-regulation; blue: significant down-regulation; black: no significant change in abundance. Brighter colors indicate stronger evidence in favor of differential abundance.

4.5 Sample size calculation for a future experiment

As in Section 3.5. Example code for the controlled spike-in experiment is given below.

```
> result.sample<-designSampleSize(data=DDAquant,numSample=TRUE,numPep=2,numTran=1,
+                               desiredFC=c(1.05,1.25),FDR=0.01,power=0.95,
+                               interference=FALSE,equalFeatureVar=FALSE)
> designSampleSizePlots(data=result.sample)
```

4.6 Quantification of protein abundance in individual samples or conditions

As in Section 3.6. Example code for the controlled spike-in experiment is given below.

```
> quantification(DDAquant)
```

5 Example: label-free DIA, a group comparison study of *S. Pyogenes*

5.1 Experimental design

This example dataset was obtained from a group comparison study of *S. Pyogenes*. Two conditions, *S. Pyogenes* with 0% and 10% of human plasma added (denoted Strep0 and Strep10), were profiled in two replicates, in the label-free mode, with a SWATH-MS-enabled AB SCIEX TripleTOF 5600 System. The identification and quantification of spectral peaks was assisted by a spectral library, and was performed using OpenSWATH software <http://proteomics.ethz.ch/openswath.html>. For reasons of space, the example dataset only contains two proteins from this study. Protein FabG shows strong evidence of differential abundance, while protein Probable RNA helicase exp9 only shows moderate evidence of differential abundance between conditions.

5.2 Reading the data

The dataset is stored as part of the package in the “long” format, as a data.frame labeled `DIARawData`. It can be accessed once the package is installed and loaded in R.

```
> head(DIARawData)
```

| | ProteinName | PeptideSequence | PrecursorCharge | FragmentIon | ProductCharge |
|-------|---------------|------------------|-----------------|-------------|---------------|
| 1 RNA | helicase exp9 | ASPIQEMTIPLALEGK | 2 | b9 | 77414 |
| 2 RNA | helicase exp9 | ASPIQEMTIPLALEGK | 2 | y10 | 77411 |
| 3 RNA | helicase exp9 | ASPIQEMTIPLALEGK | 2 | y11 | 77412 |
| 4 RNA | helicase exp9 | ASPIQEMTIPLALEGK | 2 | y14 | 77410 |
| 5 RNA | helicase exp9 | ASPIQEMTIPLALEGK | 2 | y7 | 77409 |
| 6 RNA | helicase exp9 | ASPIQEMTIPLALEGK | 2 | y9 | 77413 |

| | IsotopeLabelType | Condition | BioReplicate | Run | Intensity |
|---|------------------|-----------|--------------|-----|-----------|
| 1 | L | Strep 0% | 1 | 1 | 9747 |
| 2 | L | Strep 0% | 1 | 1 | 1272 |
| 3 | L | Strep 0% | 1 | 1 | 1295 |
| 4 | L | Strep 0% | 1 | 1 | 2332 |
| 5 | L | Strep 0% | 1 | 1 | 8178 |
| 6 | L | Strep 0% | 1 | 1 | 1737 |

5.3 Pre-processing data and quality control of MS runs

As in Section 3.3. A distinctive characteristic of DIA experiments is a potentially large number of spectral features, many of which can be reasonably expected to be non-differentially abundant between conditions. Therefore, quantile normalization is a good option for these experiments. (In this example study of *S. Pyogenes*, the intensities were normalized by quantile normalization on the full dataset with many proteins. Therefore, the code below skips the normalization step in this particular example.)

```
> DIAquant<-dataProcess(DIARawData,normalization=FALSE)
```

```
Summary of Features :
count
# of Protein      2
```

```
# of Peptides/Protein    16-25
# of Transitions/Peptide 5-6
```

Summary of Samples :

| | Strep 0% | Strep 10% |
|----------------------------|----------|-----------|
| # of MS runs | 2 | 2 |
| # of Biological Replicates | 2 | 2 |
| # of Technical Replicates | 1 | 1 |

```
> dataProcessPlots(data=DIAquant,type="ProfilePlot",ylimUp=20, ylimDown=4,
+                   featureName="NA",width=7, height=7, address=FALSE)
> dataProcessPlots(data=DIAquant,type="QCPlot",ylimUp=23, ylimDown=2,
+                   width=7, height=7, address=FALSE)
> dataProcessPlots(data=DIAquant,type="ConditionPlot",ylimUp=23, ylimDown=2,
+                   width=7, height=7, address=FALSE)
```

Figure 15 shows the profile plots for the two example proteins.




Figure 15: Profile plots for the two example proteins, Probable RNA helicase exp9 and FabG, in the study of *S. Pyogenes*. The proteins have a large number of features. The profiles show some non-parallel patterns (i.e. interferences). Features with low intensities tend to be more interference-prone.

5.4 Model-based inference

5.4.1 Setting up the model and testing proteins for differential abundance

As in Section 4.4. Since this is a label-free experiment, use `labeled=FALSE`. Example of the code and of the output with the default options for the study of *S. Pyogenes* are given below.

```
> unique(DIAquant$GROUP_ORIGINAL)

[1] Strep 0% Strep 10%
Levels: Strep 0% Strep 10%

> comparison<-matrix(c(-1,1),nrow=1)
> row.names(comparison)<-c("Strep10%-0%")
```



```
> testResultDIA<-groupComparison(contrast.matrix=comparison,data=DIAquant,labeled=FALSE)
> testResultDIA$ComparisonResult
```

| | Protein | Label | log2FC | SE | Tvalue | DF |
|---|-------------------|-------------|-------------|------------|------------|-----|
| 1 | FabG | Strep10%-0% | -1.84392847 | 0.05916153 | -31.167696 | 190 |
| 2 | RNA helicase exp9 | Strep10%-0% | -0.04755032 | 0.03005990 | -1.581852 | 296 |
| | pvalue | adj.pvalue | | | | |
| 1 | 0.0000000 | 0.0000000 | | | | |
| 2 | 0.1147512 | 0.1147512 | | | | |

5.4.2 Verifying the assumption of the model

As in Section 4.4.2. Example code for the study of *S.Pyogenes* is given below. Examples of the resulting residual plot are in Figure 16.

```
> modelBasedQCPlots(data=testResultDIA$ModelQC,type="ResidualPlots",
+                   featureName=FALSE,address=FALSE)
> modelBasedQCPlots(data=testResultDIA$ModelQC,type="QQPlots",
+                   feature.QQPlot="all",address=FALSE)
> modelBasedQCPlots(data=testResultDIA$ModelQC,type="QQPlots",
+                   feature.QQPlot="byFeature",address=FALSE)
```

For label-free DIA experiments, the modeling strategy is similar to label-free DDA experiments in Section 4.4. The code below shows the example of alternative modeling, that specifies unequal variance (`equalFeatureVar=FALSE`), and no systematic interferences (`interference=FALSE`). The change in the assumptions affected the significance status of the protein Probable RNA helicase exp9. Currently the option `equalFeatureVar=FALSE` is only available with the reduced scope of biological and technical replication.

```
> testResultDIAunequal<-groupComparison(contrast.matrix=comparison,data=DIAquant,
+                                       labeled=FALSE,equalFeatureVar=FALSE,interference=FALSE)
> testResultDIAunequal$ComparisonResult
```

| | Protein | Label | log2FC | SE | Tvalue | DF |
|---|-------------------|-------------|------------|------------|------------|-----|
| 1 | FabG | Strep10%-0% | -1.8470670 | 0.07086990 | -26.062787 | 285 |
| 2 | RNA helicase exp9 | Strep10%-0% | -0.0858471 | 0.03217217 | -2.668366 | 444 |
| | pvalue | adj.pvalue | | | | |
| 1 | 0.000000000 | 0.000000000 | | | | |
| 2 | 0.007901135 | 0.007901135 | | | | |

5.4.3 Visualizing the results of protein-level tests for differential abundance

As in Section 3.4.3. Example code for the study of *S.Pyogenes* is given below. Since at least two comparisons between conditions are needed to produce a heatmap, no heatmap is shown in this dataset. However heatmaps can be produced with other multi-condition DIA studies. Examples of the resulting comparison plots are as in Section 3.4.3.

```
> groupComparisonPlots(data=testResultDIAunequal$ComparisonResult,
+                      type="VolcanoPlot",address=FALSE)
> groupComparisonPlots(data=testResultDIAunequal$ComparisonResult,
+                      type="ComparisonPlot",address=FALSE)
```




Figure 16: (a) Residual plot for protein FabG in the study of *S. Pyogenes*, obtained under the assumption constant variance across features. The variance of one feature (colored in blue) exceeds by far the variance of the remaining features. (b) Normal quantile-quantile plot shows the pattern that is more likely due to deviations from the assumption of constant variance, and not necessarily from the assumption of Normality. (c) The diagnostics of deviations from the model assumptions per feature are typically difficult when the sample size is small.

5.5 Sample size calculation for a future experiment

As in Section 3.5. Example code for the study of *S. Pyogenes* is given below.

```
> result.sample <- designSampleSize(data=DIAquant, numSample=TRUE, numPep=5, numTran=10,
+   desiredFC=c(1.05, 1.25), FDR=0.01, power=0.95, interference=FALSE, equalFeatureVar=FALSE)
> designSampleSizePlots(data=result.sample)
```

5.6 Quantification of protein abundance in individual samples or conditions

As in Section 3.6. Example code for the study of *S. Pyogenes* is given below.

```
> quantification(DIAquant, interference=FALSE, equalFeatureVar=FALSE)
```

6 Formal description of the statistical models

6.1 SRM with stable isotope labeled reference peptides

Linear mixed effect model for SRM with stable isotope labeled reference peptides have been introduced in Chang *et al.* (2012), and are summarized in Figure 17 and Figure 18. In the notation of the figure, $i = 1, \dots, I$ is the index of a *feature*, $j = 0, \dots, J$ is the index of a *condition* or of a time point (condition 0 denotes labeled reference peptides), $k = 0, \dots, K$ is the index of a biological replicate (here called *subject*, 0 denotes labeled reference peptides), and $l = 1, \dots, L$ is the index of *run*. In presence of technical replicates, *run* simultaneously reflects the technical replication. μ_{1001} is the expected log-intensity of the first feature, first condition (i.e. with index 0), first biological replicate (i.e. with index 0) and first run, which serve as a reference. σ_{Error}^2 is the variance of the measurement error, σ_S^2 is the between-subject variance in the underlying population, $\sigma_{T \times S}^2$ is the variance due to the random between-subject interferences, σ_R^2 is the between-run variance in the underlying population, $\sigma_{F \times R}^2$ is the variance due to the random between-run interferences. A separate model is specified for each protein.

Some special cases require simpler model instances. In particular, if a protein is only represented by a single feature but has technical replicates, all the terms including *feature* are removed. If a protein is only represented by a single feature and has no technical replicates, additional terms are removed (specifically, the term *subject* for group comparison, the interaction term *time* \times *subject* for time course design, and the interaction term *condition* \times *subject* for paired design). In experiments with a single subject per condition, there is no difference between the statistical models for case-control and time-course experimental design. The terms including *subject* are removed. In absence of technical replicate the interaction term is also removed, and the simplest possible model is used.

In presence of missing values some model terms become inestimable, and the model for that protein needs to be simplified. For example, when a feature is missing completely in a condition (or a time point), the interaction *feature* \times *condition* for group comparison and paired designs (or *feature* \times *time* for time course designs) is removed. Similarly, if a feature is missing completely in a MS run for endogenous intensities, then the interaction *feature* \times *run* is removed.

If unequal variance is specified, the model is fitted using iteratively re-weighted least squares. The model terms are not changed, however the procedure affects the estimated values and the downstream statistical inference. First, fit the model with constant variance and without any weight. Next, fit a loess curve to the absolute residuals of the model above, as function of the predicted values of the model above. Calculate the weight of each peak intensity as the inverse of the fitted values from loess fit squared. Finally, re-fit the model with the new weights, and repeats the entire procedure 3 times. Currently unequal variance is only implemented with restricted scope of biological and technical replication.

| Deviation from the reference due to | | | | | | | | | | | | | | | | |
|---|---|--------------------------------------|-----------|---------------------------|---|----------------------|---------------------------------------|-------|-----------------------------------|----------------------------|--|--|--|--|--|--|
| | log(peak intensity) | Expected = reference abundance | + feature | + condition or time | + between- condition interference | + biol. replicate | + between- subject interference | + run | + between- run interference | + Random meas. error | | | | | | |
| General case | Group comparison: $y_{ijkl} = \mu_{1001} + F_i + C_j + (F \times C)_{ij} + S(C)_k + R_l + (F \times R)_{il} + \varepsilon_{ijkl}$ | | | | | | | | | | | | | | | |
| | Time course: $y_{ijkl} = \mu_{1001} + F_i + T_j + (F \times T)_{ij} + S_k + (T \times S)_{jk} + R_l + (F \times R)_{il} + \varepsilon_{ijkl}$ | | | | | | | | | | | | | | | |
| | Paired design: $y_{ijkl} = \mu_{1001} + F_i + C_j + (F \times C)_{ij} + S_k + (C \times S)_{jk} + R_l + (F \times R)_{il} + \varepsilon_{ijkl}$ | | | | | | | | | | | | | | | |
| Single feature with technical replicates | Group comparison: $y_{1jkl} = \mu_{1001} + C_j + S(C)_k + R_l + \varepsilon_{1jkl}$ | | | | | | | | | | | | | | | |
| | Time course: $y_{1jkl} = \mu_{1001} + T_j + S_k + (T \times S)_{jk} + R_l + \varepsilon_{1jkl}$ | | | | | | | | | | | | | | | |
| | Paired design: $y_{1jkl} = \mu_{1001} + C_j + S_k + (C \times S)_{jk} + R_l + \varepsilon_{1jkl}$ | | | | | | | | | | | | | | | |
| Single feature, no technical replicates | Group comparison: $y_{1jkl} = \mu_{1001} + C_j + R_l + \varepsilon_{1jkl}$ | | | | | | | | | | | | | | | |
| | Time course: $y_{1jkl} = \mu_{1001} + T_j + S_k + R_l + \varepsilon_{1jkl}$ | | | | | | | | | | | | | | | |
| | Paired design: $y_{1jkl} = \mu_{1001} + C_j + S_k + R_l + \varepsilon_{1jkl}$ | | | | | | | | | | | | | | | |
| Single subject with technical replicates | Group comparison: $y_{ij1l} = \mu_{1001} + F_i + C_j + (F \times C)_{ij} + R_l + (F \times R)_{il} + \varepsilon_{ij1l}$ | | | | | | | | | | | | | | | |
| | Time course: $y_{ij1l} = \mu_{1001} + F_i + T_j + (F \times T)_{ij} + R_l + (F \times R)_{il} + \varepsilon_{ij1l}$ | | | | | | | | | | | | | | | |
| | Paired design: $y_{ij1l} = \mu_{1001} + F_i + C_j + (F \times C)_{ij} + R_l + (F \times R)_{il} + \varepsilon_{ij1l}$ | | | | | | | | | | | | | | | |
| Single subject, no technical replicates | Group comparison: $y_{ij1l} = \mu_{1001} + F_i + C_j + R_l + \varepsilon_{ij1l}$ | | | | | | | | | | | | | | | |
| | Time course: $y_{ij1l} = \mu_{1001} + F_i + T_j + R_l + \varepsilon_{ij1l}$ | | | | | | | | | | | | | | | |
| | Paired design: $y_{ij1l} = \mu_{1001} + F_i + C_j + R_l + \varepsilon_{ij1l}$ | | | | | | | | | | | | | | | |

Figure 17: Linear mixed effect model for SRM experiments with stable isotope labeled reference peptides.

Distributional assumptions:

(a) equal variance : $\varepsilon_{ijk} \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_{Error}^2)$

(b) unequal variance (function of the expected value) : $\varepsilon_{ijk} \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_{Error}^2(\hat{y}_{ijkl}))$

Constraints:

$F_1 = 0$; $C_0 = (F \times C)_{i0} = (F \times C)_{1j} = 0$; $T_0 = (F \times T)_{i0} = (F \times T)_{1j} = 0$;

Scope of biological replication:

(a) reduced scope of biological replication: $S_0 = 0$; $(T \times S)_{j0} = (T \times S)_{0k} = 0$

(b) expanded scope of biological replication: $S_k \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_S^2)$; $(T \times S)_{jk} \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_{T \times S}^2)$

Scope of technical replication:

(a) reduced scope of technical replication: $R_1 = 0$; $(F \times R)_{i1} = (F \times R)_{1j} = 0$

(b) expanded scope of technical replication: $R_l \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_R^2)$; $(F \times R)_{il} \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_{F \times R}^2)$

Figure 18: Assumption, constraints, and scope of conclusion of linear mixed effects model for for SRM experiments with stable isotope labeled reference peptides.

6.2 Label-free DDA

Clough *et al.* (2012) introduced linear mixed effect model for label-free DDA experiments. They are summarized in Figure 19 and Figure 20. In the notation of the figure, $i = 1, \dots, I$ is the index of a *feature*, $j = 1, \dots, J$ is the index of a *condition* or of a time point, and $k = 1, \dots, K$ is the index of a biological replicate (here called *subject*). μ_{1111} is the expected log-intensity of the first feature, first condition, first biological replicate and first run, which serve as a reference. σ_{Error}^2 is the variance of the measurement error, σ_S^2 the between-subject variance in the underlying population, and $\sigma_{T \times S}^2$ the variance due to the random between-subject interferences. Experiments with technical replicates are represented by a same model (in statistical terms, the error from technical replicates is pooled into the random error). A separate model is specified for each protein.

Some special cases require simpler model instances. For some special cases, the simpler model is used. In particular, if a protein is only represented by a single feature but has technical replicates, the terms including *feature* are removed. If a protein is only represented by a single feature and has no technical replicates, additional terms are removed (specifically, the term *subject* for group comparison, the interaction term *time* \times *subject* for time course design, and the interaction term *condition* \times *subject* for paired design). In experiments with a single subject per condition, there is no difference between the statistical models for case-control and time-course experimental design. The terms including *subject* are removed. In absence of technical replicate the interaction term is also removed, and the simplest possible model is used.

As in label-based experiments, some model terms are inestimable in presence of missing values, and the model needs to be simplified for this particular protein. For example, when a feature is completely missing in one condition or one time point, the interaction *feature* \times *condition* for group comparison and paired designs (or *feature* \times *time* for time course designs) are removed.

The model with unequal variance is as in Section 6.1.

6.3 Label-free DIA

Same as in Section 6.2

| | | Deviation from the reference due to | | | | | | | |
|--|-------------------------------|-------------------------------------|-----------|---------------------|----------------------------------|------------------------|--------------------------------|------------------------|--|
| | $\log(\text{peak intensity})$ | Expected reference abundance | + feature | + condition or time | + between-condition interference | + biol. replicate | + between-subject interference | + Random meas. error | |
| General case | Group comparison: | | | | | | | | |
| | y_{ijkl} | $= \mu_{1111}$ | $+ F_i$ | $+ C_j$ | $+ (F \times C)_{ij}$ | $+ S(C)_k$ | $+ \varepsilon_{ijkl}$ | | |
| | Time course: | | | | | | | | |
| | y_{ijkl} | $= \mu_{1111}$ | $+ F_i$ | $+ T_j$ | $+ (F \times T)_{ij}$ | $+ S_k$ | $+ (T \times S)_{jk}$ | $+ \varepsilon_{ijkl}$ | |
| Single feature with technical replicates | Paired design: | | | | | | | | |
| | y_{ijkl} | $= \mu_{1111}$ | $+ F_i$ | $+ C_j$ | $+ (F \times C)_{ij}$ | $+ S_k$ | $+ (C \times S)_{jk}$ | $+ \varepsilon_{ijkl}$ | |
| | | | | | | | | | |
| Single feature, no technical replicates | Group comparison: | | | | | | | | |
| | y_{1jkl} | $= \mu_{1111}$ | $+ C_j$ | $+ S(C)_k$ | $+ \varepsilon_{1jkl}$ | | | | |
| | Time course: | | | | | | | | |
| | y_{1jkl} | $= \mu_{1111}$ | $+ T_j$ | $+ S_k$ | $+ (T \times S)_{jk}$ | $+ \varepsilon_{1jkl}$ | | | |
| Single feature, no technical replicates | Paired design: | | | | | | | | |
| | y_{1jkl} | $= \mu_{1111}$ | $+ C_j$ | $+ S_k$ | $+ (C \times S)_{jk}$ | $+ \varepsilon_{1jkl}$ | | | |
| | | | | | | | | | |
| Single subject with technical replicates | Group comparison: | | | | | | | | |
| | y_{ij1l} | $= \mu_{1111}$ | $+ F_i$ | $+ C_j$ | $+ (F \times C)_{ij}$ | $+ \varepsilon_{ij1l}$ | | | |
| | Time course: | | | | | | | | |
| | y_{ij1l} | $= \mu_{1111}$ | $+ F_i$ | $+ T_j$ | $+ (F \times T)_{ij}$ | $+ \varepsilon_{ij1l}$ | | | |
| Single subject, no technical replicates | Paired design: | | | | | | | | |
| | y_{ij1l} | $= \mu_{1111}$ | $+ F_i$ | $+ C_j$ | $+ (F \times C)_{ij}$ | $+ \varepsilon_{ij1l}$ | | | |
| | | | | | | | | | |
| Single subject, no technical replicates | Group comparison: | | | | | | | | |
| | y_{ij11} | $= \mu_{1111}$ | $+ F_i$ | $+ C_j$ | $+ \varepsilon_{ij11}$ | | | | |
| | Time course: | | | | | | | | |
| | y_{ij11} | $= \mu_{1111}$ | $+ F_i$ | $+ T_j$ | $+ \varepsilon_{ij11}$ | | | | |
| Single subject, no technical replicates | Paired design: | | | | | | | | |
| | y_{ij11} | $= \mu_{1111}$ | $+ F_i$ | $+ C_j$ | $+ \varepsilon_{ij11}$ | | | | |
| | | | | | | | | | |

Figure 19: Linear mixed effect model for label-free DDA experiments.

Distributional assumptions:

(a) constant variance : $\varepsilon_{ijk} \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_{Error}^2)$

(b) unequal variance (feature-dependent variance) : $\varepsilon_{ijk} \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_{Error}^2(\hat{y}_{ijkl}))$

Constraints:

$F_1 = 0$; $C_1 = (F \times C)_{i1} = (F \times C)_{1j} = 0$; $T_1 = (F \times T)_{i1} = (F \times T)_{1j} = 0$;

Scope of biological replication:

(a) reduced scope of biological replication: $S_1 = 0$; $(T \times S)_{j1} = (T \times S)_{1k} = 0$

(b) expanded scope of biological replication: $S_k \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_S^2)$; $(T \times S)_{jk} \stackrel{iid}{\sim} \mathcal{N}(0, \sigma_{T \times S}^2)$

Figure 20: Assumption, constraints, scope of conclusion of linear mixed effects model for label-free DDA experiments.

References

- Amaratunga, D. and Cabrera, J. (2001). Analysis of Data From Viral DNA Microchips. *Journal of the American Statistical Association*, **96**(456), 1161–1170. [14](#), [35](#)
- Benjamini, Y. and Hochberg, Y. (1955). Controlling the false discovery rate: a practical and powerful approach to multiple testing. *J.R. Statist. Soc. B*, **57**(1), 289–300. [21](#)
- Chang, C., Picotti, P., Hüttenhain, R., Heinzelmann-Schwarz, V., Jovanovic, M., Aebersold, R., and Vitek, O. (2012). Protein significance analysis in selected reaction monitoring (SRM) measurements. *Mol. and Cell. Proteomics*, **11**, 1–12. [4](#), [43](#)
- Clough, T., Thaminy, S., Ragg, S., Aebersold, R., and Vitek, O. (2012). Statistical protein quantification and significance analysis in label-free LC-MS experiments with complex designs. *BMC Bioinformatics*, **13**(Suppl 16), S6. [4](#), [46](#)
- Cox, J. and Mann, M. (2008). MaxQuant enables high peptide identification rates, individualized p.p.b.-range mass accuracies and proteome-wide protein quantification. *Nature Biotechnology*, **26**(12), 1367–1372. [3](#)
- Gatto, L. and Lilley, K. S. (2012). MSnbase-an R/Bioconductor package for isobaric tagged mass spectrometry data visualization, processing and quantitation. *Bioinformatics*, **28**, 288–289. [3](#)
- Kutner, M., Nachtsheim, C., Neter, J., and Li, W. (2005). *Applied linear statistical models*. McGraw Hill, New York. [19](#)
- MacLean, B., Tomazela, D. M., Shulman, N., Chambers, M., Finney, G., Frewen, B., Kern, R., Tabb, D. L., Liebler, D. C., and MacCoss, M. J. (2010). Skyline: An open source document editor for creating and analyzing targeted proteomics experiments. *Bioinformatics*, **26–27**, 966. [4](#)
- Mueller, L. N., Rinner, O., Schmidt, A., Letarte, S., Bodenmiller, B., Brusniak, M., Vitek, O., Aebersold, R., and Müller, M. (2007). SuperHirn - a novel tool for high resolution LC-MS-based peptide/protein profiling. *Proteomics*, **7**, 3470–3480. [3](#), [34](#)
- Picotti, P., Bodenmiller, B., Mueller, L. N., Domon, B., and Aebersold, R. (2009). Full dynamic range proteome analysis of *S. cerevisiae* by targeted proteomics. *Cell*, **138**, 795–806. [13](#)
- Sturm, M., Bertsch, A., Gröpl, C., Hildebrandt, A., Hussong, R., Lange, E., Pfeifer, N., Schulz-Trieglaff, O., Zerck, A., Reinert, K., and Kohlbacher, O. (2008). OpenMS – An open-source software framework for mass spectrometry. *BMC Bioinformatics*, **9**(1), 163. [3](#)
- Surinova, S., Hüttenhain, R., Chang, C., Espona, L., Vitek, O., and Aebersold, R. (2013). Automated selected reaction monitoring data analysis workflow for large-scale targeted proteomic studies. *Nature Protocols*, **8**, 1602–1619. [4](#)